**HEMOSAL****GB****Description**

Astringent gel - 25% aluminum sulfate

**Avantages**

- Available in pre-filled syringe: hygienic and ergonomic.
- Makes initial cord packing easier by providing lubrication: allows the cord to glide into the sulcus.
- No tissue necrosis or blackening of tissue. compatible with all impression materials.
- Reduced tissue trauma.
- Blue in colour for easy visibility and placement.
- Plesant raspberry flavour.
- The viscous gel, unlike runny liquids or syrups, remains in place for the extent of the procedure.
- Versatile in uses.

**Indications**

- To help in gingival retraction, with the appropriate retraction cord.
- Help in control of bleeding during impression procedures.
- Help in reducing gingival inflammation.
- Great of class V restorations or if tissue is cut during composite procedures.
- Other tissues management.
- In endodontic procedures, when a canal is bleeding.

**Instruction for use**

- Put a disposable applicator on the syringe.
- Impression
- Gently apply a thin layer of gel directly to the sulcus area.
- Select the appropriate cord (00, 0 or 01).
- Pack the selected cord through the gel, into the sulcus.
- Waite 2 to 6 minutes.
- Remove the cord, rinse and dry.
- Take the impression as usual.

**Parodontology & dentistry**

- Apply Hemosal with pressure, directly on haemorrhaging area, with a bristle brush or a Elsobrush.
- Wait for 2 to 4 minutes and rinse.

**Endodontic treatments**

- Apply a paper cone into the gel and apply the prepared cone in the haemorrhagic canal. For 20 to 30 seconds. Reapply another paper cone if necessary.

**Contraindications**

The use of Hemosal is contraindicated if the patient is known to be allergic to any of the ingredients in Hemosal.

**Important notice**

- Do not throw away the box and the instruction for use until the product is completely finished.
- High humidity could damage the packaging on which important instructions and legal information are printed.
- Storage at a too high temperature could cause to age prematurely and inhibit its polymerisation.
- Storage at a too low temperature leads to the decomposition of the product which will render it completely unusable.

**Disclaimer**

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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 Conserver au sec  
 Keep in dry place

 Température de stockage  
 Temperature of storage

 Attention, voir notice d'instructions  
 Prior to use, carefully read instructions of use

 Adresse du fabricant  
 Address of manufacturer

# Hemosal

## ELSODENT

**G-PHARMA**  
 35, avenue du Gros Chêne  
 B.P. 10279 HERBLAY  
 95617 CERGY-PONTOISE  
 Cedex FRANCE  
 elsodent@elsodent.com  
 www.elsodent.com