PureFill 2

GB

Definition

Light-cure nano-hybrid composite, for anterior and posterior restorations. Highly filled. Restorations of class I, II, III, IV and V.

Safety notes

Warnings

- PureFill 2 contains polymerizable methacrylate resins which may irritate skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
- Avoid eye contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- Avoid skin contact to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.
- Avoid contact with oral soft tissues/mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosa persists, seek medical attention.

Precautions

- Wear suitable protective eyewear, mask, clothing and gloves.
- Protective eyewear is recommended for patients.

Contraindications

PureFill 2 is contraindicated in patients with a history of allergic reaction to methacrylate resins or any other of the components. Do not use PureFill 2 in cases of pulpar exposition, pulpitis.

Instruction for use

1. Select the appropriate VITA shade with day-light.
2. Prepare the tooth cavity in a conventional manner. The margins are prepared by slightly beveling the occlusal surfaces and producing slight flares on interproximal cavities.
3. For classes II, III and IV, a matrix band is placed and secured by wedges; it should be loosened to allow full contours when packing the restorative material into the cavity.
4. If the cavity is closed to the pulp, apply a eugenol-free cavity liner (ex: calcium hydroxide, antibacterial liner).
5. Etch the enamel and dentin.
6. Apply the light curing bonding agent of your choice, following manufacturer's instructions.
7. Apply PureFill 2 and, using a high-power light curing unit (>1200 mW/cm²), polymerize PureFill 2 in 3 mm layer. Exposure time is 15 seconds for light shades (A1, A2, C2, B2 etc) and 20 seconds for dark shades (A3, A3.5, A4, etc).
8. Finish and polish the restoration.

Storage & shelf-life

Do not expose the product to sun-light or warm sources. Store at a temperature not exceeding 25°C. When stored under such conditions, PureFill 2 has a shelf-life of 36 months.

Presentations & references

PUR2-3* - 3g syringe. Available in the following VITA shades: A1/B1, A2, A3, A3.5, A4, B2, B3, C2 C3, C4, D3, Incisal.

PUR2-5* - Box of 20 x 0.25 g compules, available in the following VITA shades : A1/B1, A2, A3, A3.5, B2. *mentioned the shade.

Important notice

- Do not throw away the box and the instruction for use until the product is completely finished.
- High humidity could damage the packaging on which important instructions and legal information are printed.
- Storage at a too high temperature could cause to age prematurely and inhibit its polymerization.
- Storage at a too low temperature leads to the decomposition of the product which will render it completely unusable.

Disclaimer

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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