**Kerr™**

# Maxcem Elite™ Chroma

Self-Etch / Self-Adhesive Resin Cement with Cleanup Indicator

Instructions for use

## **Maxcem Elite Chroma**

### **Description**

Maxcem Elite Chroma is a unique resin cement that has a gel-state color indicator which visually displays the optimal time to remove the excess cement, making it one of the smartest cements available on the market today. Cleaning the cement around the margins at the optimal time facilitates clean-up, ensures success of the restoration and aids in preventing post-operative complications associated with residual cement. The cement turns pink when dispensed from the syringe and the pink color will disappear once the cement reaches its gel state, indicating ideal timing for excess cement removal.

Maxcem Elite Chroma is radiopaque for easy identification in x-rays and can be used in self-cure or light cure mode, allowing for maximum flexibility in its use. Maxcem Elite Chroma is packaged in a dual barrel syringe with single-use automix tips and optional curved dispensing tips to allow user to deliver desired volume of cement directly into

restoration and/or tooth preparation.

### **Indications for Use**

Maxcem Elite Chroma is indicated for cementation of all indirect restorations including ceramic, resin, and metal-based inlays/onlays, crowns, bridges, posts, and veneers\*. Additional indications include cementation of crowns to implants.

\*Adhesive application on the preparation is required for veneer cementation using Maxcem Elite Chroma. Please see the "Veneer Restorations" section for more details.

### **Contraindications**

1. This product shall not be used on patients who are allergic to methacrylate resin as it may cause contact dermatitis.
2. This product shall not be placed in contact with pulp tissue. Direct exposure of this product may harm the pulp.

### **Warnings**

1. Uncured methacrylate resin may cause contact dermatitis and damage the pulp. Avoid contact with skin, eyes, and soft tissue. Wash thoroughly with water after contact.

### **Precautions**

1. Maintain restoration with moderate pressure while cleaning excess cement.
2. Allow the cement to fully cure before finishing and polishing or adjusting occlusion to avoid potential debonding.
3. Do not use phosphoric acid to etch dentin when using Maxcem Elite Chroma, as this may cause post-operative sensitivity.
4. Do not treat or clean the zirconia-based or metal-based restorations with phosphoric acid or hydrofluoric acid (HF) which can negatively impact adhesion.
5. Do not use hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) as its residues can adversely affect the setting reaction and adhesion of Maxcem Elite Chroma.

6. Do not use provisional cements containing eugenol, as these materials can inhibit the polymerization of Maxcem Elite Chroma.
7. Use a rubber dam to keep the preparation or cavity free of any oral fluid contamination during cementation.
8. Do not open foil package until initial use. Use the product within 6 months of opening the pouch but always ensure that it is not past the expiration date printed on the syringe label.
9. If prolonged contact with oral tissue occurs, flush with copious amount of water.

### **Adverse reactions**

None known.

### **Instructions for Use**

**Note: To prevent cross-contamination between patients, automix tips are for single patient use only.**

#### **Preparation of the Restoration Surface**

Prepare the internal surface of the restoration according to manufacturer's instructions. Following are some general recommendations:

**Etchable Ceramic/Porcelain Restorations:** Sandblast the internal surface with 50µm alumina at a pressure of 30 psi (0.2 MPa) and etch the restoration with hydrofluoric acid. Application of a silane primer is optional.

**Metal-based, Zirconia-based, and Composite Restorations:** Sandblast the internal surface with 50µm alumina at a pressure of 60 psi (0.4 MPa) for metal-based or zirconia-based restorations, or a pressure of about 15 psi (0.1 MPa) for composite restorations.

**Resin bonded prostheses (Maryland bridges):** Etching of enamel surface with phosphoric acid gel etchant is recommended prior to cementation with Maxcem Elite Chroma. Treat the bonding surface of the restoration according to manufacturer's recommendations.

## **Cementation of Inlay, Onlay, Post, Crowns and Bridges** **(For implant restorations and veneers, refer to the appropriate section)**

### **Preparation**

1. Remove the provisional restoration.
2. Use the OptiClean rotary instrument or your choice of cleaning method to remove provisional cement and residue from the preparation. Rinse thoroughly.
3. Air dry the preparation.
4. Seat the restoration to ensure proper fit. Adjust if needed and thoroughly clean the intaglio surfaces of the restoration after adjustment.
5. Select the desired shade of Maxcem Elite Chroma.
6. Remove the syringe from the foil pouch and record the date on the syringe label. Remove the syringe cap. Before the first use, bleed the syringe by squeezing out a small quantity of the pastes onto a pad to equalize the base and catalyst pastes in the syringe. Secure the appropriate automix tip on the dual syringe cartridge by turning the tip clockwise. The wide 4:1 automix tips are intended for use with the intraoral tips and root-canal tips to dispense Maxcem Elite Chroma. Leave the used mixing tip on the syringe as a cap until the next

application.

## Cementation

7. Dispensing of the material is done as follows:

- **Inlay/Onlay** - Dispense the cement directly into the preparation covering all surfaces. Seat the restoration with moderate pressure.
- **Crown** - Fill the crown with cement. Seat the restoration with moderate pressure.
- **Post** - Prepare the post space. Size and fit the post. Dispense the cement onto the post and directly into the canal space using the root-canal tips. Seat and rotate the post slightly upon insertion to avoid trapping air. Allow the cement to slowly flow from canal space.

8. Maxcem Elite Chroma has a gel-state color indicator which visually displays the optimal time to remove excess cement. The cement turns pink when dispensed from the syringe and the pink color will disappear once the cement reaches its gel state, indicating ideal timing for excess cement removal with a scaler or explorer. For the

best cleanup experience, follow one of the following options:

- **Option 1 – Self Cure:** Wait until the pink color in the cement fades away before initiating clean-up. This will take approximately 2 minutes, though the color change is the best indication of the gel-state. Maintain pressure on the restoration while removing excess cement, and until the cement has fully set.
  - **Option 2 – Tack Cure:** You can initiate the gel-state by tack-curing excess cement with a curing light for approximately 2-3 seconds (pink color will disappear immediately). Excess cement should be immediately removed with a scaler or explorer. Maintain pressure on the restoration while removing excess cement.
  - **Option 3 – Immediate Removal:** Remove all uncured pink-colored excess cement using a sponge pellet or microbrush while maintaining pressure on the restoration. Cover the margins with Kerr NX3 Try-in gel or a glycerin gel of your choice. Immediately light-cure all surfaces including margins using Demi Ultra curing light for 10 seconds or allow the cement to self-cure for at least 4 minutes after seating of the restoration. Rinse glycerin gel off with water. Proceed directly to step 10.
9. After removal of the excess cement, light cure all surfaces including margins using Demi Ultra curing light for 10

seconds. See footnote for recommended curing times for other curing light models.\*\*

10. For metal-based and opaque restorations such as PFM and Zirconia, it will not be possible to light cure the cement. Before proceeding to finishing and polishing, ensure that the cement has been allowed to cure for at least 4 minutes after seating.
11. If needed, use finishing strips for proximal areas. Polish margins with Gloss Plus discs, cups, or points. Use Hi-Luster Plus polishers for a final high luster shine.
12. Verify the occlusion.

For additional guidance, refer to the printed technique guide included in the packaging.

## Implant Restorations

1. Prepare the implant abutment per manufacturer's instructions.
2. If needed, clean the abutment and air dry.
3. Seat the restoration to ensure proper fit. After try-in and occlusal adjustment, thoroughly clean the bonding

surface of the restoration.

4. Select the desired shade of Maxcem Elite Chroma.
5. **Follow the steps 7 – 12 in the “Cementation” section for crowns above.**

## Veneer Restorations

For cementation of veneers (ideal for 1-2 units at a time), an adhesive must first be applied to the tooth. For best results, select OptiBond™ Solo Plus† Adhesive or OptiBond XTR†† Self-Etch Universal Adhesive. Consult adhesive product instructions for use.

After adhesive is applied and cured, follow the instructions for cementation:

1. Select the desired shade of cement.
2. Remove the syringe cap. Before the first use, bleed the syringe by squeezing out a small quantity of pastes from each side onto a pad to equalize the base and catalyst paste in the syringe. Secure the appropriate automix tip on the dual syringe cartridge by turning the tip clockwise.

3. Dispense the cement directly onto the internal surface of the veneer. Place the veneer gently onto the tooth, allowing the cement to slowly flow from all margins.
4. Spot cure the veneer into place on the facial surface, away from the margins, using a small diameter light guide for 10 seconds.
5. After removing the excess cement, light cure all surfaces using a Demi Ultra curing light for 10 seconds.\*\*

**Additional Notes**

†for OptiBond S

††for OptiBond Versa

**Note: Maxcem Elite Chroma minimum working time is 1 minute 20 seconds at room temperature; Set time is less than 4 minutes at oral temperature. These times may vary based on storage conditions, temperature, humidity, and age of the product.**

**\*\* Recommended Curing Times:**

**Demi Ultra, Demi Plus/Demi, and L.E.Demetron II: 10 seconds.**

**L.E.Demetron I & Optilux 501: 20 seconds.**

**For all other lights, see manufacturer's recommendations.**

**STORAGE:** Store the product within a temperature range of 4°C (39°F) to 25°C (77°F). Avoid excess heat. For optimum freshness, store product in a refrigerator. Allow the product to reach room temperature prior to use. Keep product in foil package until initial use. Once the foil pouch is opened, use within 6 months and prior to expiration date printed on syringe label.

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**Kerr Corporation**
**Kerr Italia, S.r.l.**