**EN** **FR** **DE** **IT** **ES** **TR**

# Gel Desensitizer

Desensitizer for Hypersensitive Teeth

## EN

**INSTRUCTIONS**

**IMPORTANT :** READ ALL INSTRUCTIONS THOROUGHLY BEFORE USE. KEEP THIS LEAFLET AND REFER TO IT PERIODICALLY.

Desensitizer for Hypersensitive Teeth

## Gel Desensitizer

**FOR DENTIST USE ONLY**

**Contraindications:**

Patients who have a history of hypersensitivity to this product.

**1. What is Gel Desensitizer?**

Gel Desensitizer is a single-syringe dentin tubule sealant with fluoride. It reacts chemically with tooth structure and forms a precipitate that blocks dental tubules, to prevent thermal, mechanical and chemical stimulation of odontoblastic processes. Application of Gel Desensitizer is indicated for topical treatment of dental hypersensitivity due to patent dentinal tubules, and will immediately reduce it. Treatment will remain effective for a significant time period. The use of Gel Desensitizer will NOT relieve the discomfort of pulpitis.

**Indications:**

Gel Desensitizer is used for topical treatment of hypersensitivity dentin.

### 2. Contents

Gel Desensitizer

① Gel Desensitizer ----- 3mL ----- Class IIa

② Plastic Needle Tip ----- 10 (with grey Needle Cap) --- Class I

③ Instructions

**3. Composition**

Copolymer with sulfonic acid group, Oxalic acid, Water, Sodium fluoride, pH buffer (Potassium salt), Thickener, etc.

**4. Precautions**

**4.1 Safety**

This product contains substances that may cause allergic reactions. Read the following points below.

**1) Precautions for use**

- ① Avoid mixing with other products.
- ② If accidentally swallowed during use, immediately have the patient gargle with a water, consult a physician if necessary.
- ③ Plastic Needle Tip is one-time use only.
- ④ If you push the plunger strongly, the gel might splatter around. Please push the plunger gently.

**2) Important precautions**

- ① Discontinue the use immediately if signs of irritation (rash, etc.) appear. See a physician.
- ② If an operator has a history of hypersensitivity to this product, use protective gears to avoid direct contact to it. If the hypersensitivity reactions appear, see a physician.
- ③ This product is acidic. Avoid contact with oral mucosa, skin or eyes. Contaminated skin should be thoroughly rinsed with running water. Materials left on the tissue may cause, for example, swelling of the lips or white sore at the corner of the mouth. In case of the eye contact, flush immediately with copious amounts of water and consult an ophthalmologist if necessary.
- ④ Be cautious on a risk of causing gingival inflammation by rubbing the gingiva excessively with cotton pledget and rubber cup etc.

### 4.2 To get the best results with Gel Desensitizer

**1) Storage conditions**

- ① Store Gel Desensitizer in dark at room temperature (1°C - 30°C). Avoid exposure to direct sunlight and extreme temperatures (above 30°C or below 1°C).
- ② Do not use this product after the expiration date. See outer package for expiration date.

**2) Precautions about how to use**

- ① Do not use Gel Desensitizer after the root planing if reattachment is expected.
- ② Desensitizing efficacy may be reduced if a patient has bruxism or a tooth with a large non-carious cervical lesion (wedge-shaped defect).
- ③ Once a tooth is sealed with Gel Desensitizer, the bond strength of subsequently applied adhesives will be compromised. To seal a surface that will be bonded later, use a bonding agent such as Hybrid Bond (from SUN MEDICAL Co., Ltd.) instead of Gel Desensitizer. If a surface treated with Gel Desensitizer must be bonded, expose the fresh surface by sanding or etching, Application of Green Activator or Dentin Etchant Gel (from SUN MEDICAL Co., Ltd.), rubbed for 10 seconds with a saturated pledget before bonding with Super-Bond C&B (from SUN MEDICAL Co., Ltd.), will attain good bond strengths.
- ④ This product contains fluoride. Avoid contact with the glass product as it can react / etc.
- ⑤ Recap the syringe immediately after each use. Dried out this product can clog the syringe if left uncapped.
- ⑥ Material's taste may change slightly depends on the storage environment and duration. However, it does not affect the performance.

**3) Product life time**

Product life time can be shortened depending on the circumstances of use. Read instructions of relevant products carefully before use.

### 4.3 Disposal

Dispose of empty container or package according to the law about the waste.

***As in any dental treatment, the patient's individual constitution and the unique requirement of the clinical case at hand must be considered before selecting materials and conditions for use.***

**5. How to use Gel Desensitizer**

**① Surface preparation**

Tooth surface preparation: Remove plaque and clean the tooth surface with a cotton pledget and floss thoroughly.

**② Application of Gel Desensitizer**

- a) Remove the syringe cap and attach a Plastic Needle Tip to the syringe.
- b) Apply a thin layer to the lesion or the tooth surface around it and leave it for 30 sec. Recommended quantity per tooth is the size of a grain of rice (Φ: 3mm).
- c) Attach Plastic Needle Tip to Gel Desensitizer and twist firmly until it is tight.

**③ Rinse**

Have the patient rinse the mouth with water thoroughly.

**If the desensitizing is not satisfactory with the first application, repeat the treatment at successive appointments to increase desensitization.**

\* To obtain higher desensitizing results or for the case establishing the clean surface is difficult due to the severe sensitivity, apply Gel Desensitizer with **rubbing/pumping motion for more than 5 seconds** using cotton pledget or rubber cup by following the instruction below.

- 1) Selection and preparation of applicators
Select a cotton pledget according to the size of the affected area or attach the rubber cup etc. to the contra-angle.
- 2) Rubbing/pumping motion
Dispense the desensitizer onto the lesion and/or surrounding tooth surface (standard quantity per tooth is the size of a grain of rice). Rub and pump for more than 5 seconds using the cotton pledget. If using the rubber cup, pick up the desensitizer (about the size of a grain of rice per tooth) in the cup, apply by rubbing and pumping for more than 5 seconds at ~1,000 rpm so that Gel Desensitizer does not scatter. The use of cotton pledget and rubber cup etc. is preferable for the lesions close to gingiva or interproximal areas.

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