mectron

medical technology

EN

FR

ES

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# → TABLE OF CONTENTS

EN

01 →	INTRODUCTION .....	4
01.1	Intended use of the multipiezo .....	4
01.2	Description of the device .....	5
01.3	Responsibility waiver .....	5
01.4	Safety precautions .....	5
01.5	Symbols .....	7
02 →	IDENTIFICATION DATA .....	8
02.1	Identification plate of the device .....	8
02.2	Identification data of the handpiece .....	8
02.3	Identification data of the inserts .....	8
03 →	DELIVERY .....	9
03.1	List of the components of the multipiezo .....	9
04 →	INSTALLATION .....	10
04.1	First installation .....	10
04.2	Safety precautions during installation .....	10
04.3	Connection of the accessories .....	10
05 →	USE .....	12
05.1	Switching the device on and off .....	12
05.2	Description of the keyboard .....	13
05.3	Safety precautions before and during use .....	16
05.4	Instructions for use .....	18
05.5	Important information on the inserts .....	20
06 →	FLUSH FUNCTION .....	21

- 07 → DISASSEMBLY OF PARTS FOR CLEANING AND STERILIZATION..... 22
- 08 → CLEANING AND STERILIZATION ..... 24
  - 08.1 Cleaning of the foot pedal ..... 24
  - 08.2 Cleaning of the device body ..... 25
  - 08.3 Cleaning and sterilization of the handpiece ..... 26
  - 08.4 Cleaning and sterilization of the inserts ..... 29
  - 08.5 Cleaning and sterilization of the torque wrench ..... 32
  - 08.6 Cleaning of the bottle ..... 34
- 09 → MAINTENANCE ..... 36
  - 09.1 Replacement of the peristaltic pump ..... 36
- 10 → DISPOSAL MODES AND PRECAUTIONS ..... 37
- 11 → TECHNICAL DATA ..... 38
  - 11.1 Electromagnetic compatibility EN 60601-1-2 ..... 39
- 12 → TROUBLESHOOTING..... 43
  - 12.1 Diagnostic system and symbols on the keyboard ..... 43
  - 12.2 Quick solution to problems ..... 44
  - 12.3 Replacement of the fuses ..... 46
  - 12.4 Sending the device to an Authorized Mectron Service Service .... 47
- 13 → WARRANTY ..... 48

## 01 → INTRODUCTION

EN

Carefully read this manual before proceeding with the installation, use, maintenance, or other operations on the device. Always keep this manual within reach.

**IMPORTANT:** to prevent harm to persons or damages to objects, read all the "Safety precautions" present in the manual with special care. Depending on their degree of seriousness, the safety precautions are classified with the following indications:

-  **WARNING**  
(always refers to personal injury)
-  **CAUTION**  
(refers to possible damage to property)

The purpose of this manual is to make the operator knowledgeable of the safety

precautions, the installation procedures, and the instructions for a correct use and maintenance of the device and its accessories.

Use of this manual for purposes other than those strictly tied to the installation, use and maintenance of the device is forbidden. The information and illustrations in this manual are updated as of the date of issue reported on the last page.

Mectron is committed to continuously update its products with possible modifications to device components. In case you uncover discrepancies between what is described in this manual and the device in your possession, contact your Retailer or the After-Sales Service of Mectron for clarifications and support.

### 01.1 → INTENDED USE OF THE MULTIPIEZO

By means of the appropriate inserts, the following dental applications can be performed:

- scaling: all the procedures for the removal of bacterial plaque deposits and supra-gingival, sub-gingival, and interdental calculi, and for the removal of stains;
- periodontology: periodontal treatment for the scaling and root-planing/debridement, including cleaning and irrigation of the periodontal pocket;
- cleaning of the implant site surfaces;
- endodontics: all treatments for root canals preparation, irrigation, filling, condensation of the gutta-percha, endodontic revision and retrograde preparation;
- restoration and prosthesis: cavity preparation and removal of carious tissues, removal of denture and of excess restoration material, condensation of the amalgam, finishing of the cervical margin.

 **WARNING:** The device must be used in a dental practice or ambulatory. The device cannot operate in environments where anesthetic or flammable mixtures are present.

 **WARNING: Qualified and specialized personnel.** The device must be used exclusively by specialized and appropriately trained personnel. The use of the device

does not cause side effects if it is used correctly. An improper use might cause tissues heating.

 **WARNING:** Use the device only for the intended use. Failure to observe this precaution may cause serious injuries to the patient, the operator, and damages/breakdowns to the device.

## 01.2 → DESCRIPTION OF THE DEVICE

The multipiezo is a multifunctional ultrasonic piezoelectric scaler. It has been designed in order to offer a product with innovative design and exclusive technical characteristics to the operator, and maximum comfort during the treatment to the patient.

The user interface has been optimized making all the functions readily available by

integrating them in the touch keyboard. The handpiece has a LED light 360° adjustable and can be autoclaved at 135 °C. The device has an automatic sync circuit which optimizes frequency and power for each available insert, so that the best performance is always assured.

## 01.3 → RESPONSIBILITY WAIVER

The manufacturer Mectron disclaim any liability, expressed or implied, and shall have no responsibility for any direct, indirect or other damages and personal injury arising out in connection with any errors in the use of the device and its accessories.

The manufacturer Mectron shall be under no liability, expressed or implied, with respect to any damages (personal injury and/or damage to property) which might arise or be caused, whether by the customer or by any of the users of the product and its accessories, as result of:

- 1 Use or procedures different than those specified in the intended use of the product;
- 2 The environmental conditions for the preservation and storage of the device are not compliant to the precautions indicated in the Chapter 11 - TECHNICAL DATA;
- 3 The device is not used in compliance with all the instructions and precautions described in this manual;
- 4 The electrical system in the premises in which the device is used is not compliant to the norms in force and to the relative precautions;
- 5 The assembly operations, extensions, adjustments, updates, and repairs on the device are performed by personnel not authorized by Mectron;
- 6 Improper use, mistreatments, and/or incorrect interventions;
- 7 Any and all attempts to tamper with or modify the device, under any circumstance;
- 8 Use of non-original Mectron inserts that entail a finite damage to the threading of the handpiece, thus compromising correct operation and causing risk of harm to the patient;
- 9 Use of non-original Mectron inserts, used in accordance to designed and tested settings of Mectron original inserts. The correct use of the settings is guaranteed only with original Mectron inserts;
- 10 Lack of stock materials (handpiece, inserts, wrenches) to be used in the event of device stop due to fault or of inconveniences.

## 01.4 → SAFETY PRECAUTIONS

**⚠ WARNING: Contraindications.** Do not use the multipiezo on patients who carry heart stimulators (Pace-makers) or other implantable electronic devices. This precaution also applies to the operator.

**ⓘ CAUTION: Contraindications. Ultrasonic scaler.** Do not perform treatments on prosthetic artifacts made of metal or ceramics. The ultrasonic vibrations could lead to the de-cementing of the artifacts.

**⚠ WARNING: Contraindications.** Do not

perform scaling treatments without water spray in order to avoid the insert overheating which may cause damages to the tooth. Treatments with no water spray can be carried out only with "Dry Work" inserts which do not have the water passage.

**⚠ WARNING: Contraindications. Interference from other equipment.** An electrical scalpel or other electro-surgical units near the multipiezo device may interfere with its correct operation.

**⚠ WARNING: Contraindications.** **Interference with other equipment.** Though compliant to the standard IEC 60601-1-2, the multipiezo may nonetheless interfere with other devices nearby. The multipiezo must not be used near to or stacked on other devices. However, if this were to prove necessary, you must check and monitor correct operation of the device in that configuration.

**⚠ WARNING: Do not install the device in places where there is a risk of explosion.** The device cannot operate in environments where anesthetic or flammable mixtures are present.

**⚠ CAUTION:** In the case that the end user, when operating in his or her own medical study or clinic, must subject the electro-medical equipment and systems to periodical inspections in order to adhere to imposed requirements, the test procedures that must be applied to electro-medical equipment and systems to evaluate safety must be performed in line with norm EN 62353 'Electro-medical devices - Periodical inspections and tests to be performed after repair interventions on electro-medical devices'.

**⚠ WARNING: Checking device status before the treatment.** Always check that there is no water underneath the device. Before every treatment, always check that the device works perfectly and that the accessories are efficient. In case you uncover operating abnormalities, do not perform the treatment. Contact an Authorized Mectron Service Center if the abnormalities concern the device.

**⚠ CAUTION:** The electrical system of the premises in which the device is installed and used must be compliant to the norms in force and to the relative electrical safety precautions.

**⚠ CAUTION:** To avoid any risk of electric shock this device must be grounded.

**⚠ WARNING: All the accessories of the new or repaired devices are not sterile.** All new and repaired accessories are supplied in non-sterile conditions. Before use, and after each treatment, they must be cleaned

and sterilised in strict compliance with the instructions given in the Chapter 08 - CLEANING AND STERILIZATION.

**⚠ WARNING: Infections control.** To ensure maximum safety of the patient and the operator, before using all the reusable parts and accessories, make sure you have first cleaned and sterilized them, following the instructions of the Chapter 08 - CLEANING AND STERILIZATION.

**⚠ CAUTION: Contraindications.** Allow reusable, autoclavable items (the handpiece, the inserts, the torque wrench, and any other accessory that can be sterilized ) to gradually return to room temperature after steam sterilization and prior to usage. The cooling process must not be accelerated.

**⚠ WARNING: Breakage and wear-out of the inserts.** High frequency oscillations and wear-out may, in rare circumstances, lead to the breakage of the insert.

Deformed or otherwise damaged inserts are susceptible to breakage during their use. These inserts must never be used. If an insert breaks, check that none of its fragments remain in the treated part and, at the same time, apply effective suction to remove them.

The patient must be instructed to breathe through his nose during the treatment, or a dental dam must be used to prevent the patient from ingesting fragments of broken inserts.

When the nitride coating wears out, the cutting efficiency decreases; re-sharpening the insert damages it and is therefore forbidden. Check that the insert is not worn out.


During the intervention, frequently check that the insert is intact, especially in its apical part.

During the intervention, avoid prolonged contact with retractors or with metallic instrumentation in use. Do not exert excessive pressure on the inserts during their use.

**⚠ WARNING:** Only use original Mectron inserts, accessories, and spare parts.

**⚠ CAUTION:** No modification of this equipment is allowed.

## 01.5 → SYMBOLS

	Serial number
	Lot number
	Product code
	CAUTION: read the instructions for use
	Operating instructions
	Temperature limitation - transport and storage conditions
	Humidity limitation - transport and storage conditions
	Atmospheric pressure limitation - transport and storage conditions
	It indicates compliance with the CE 93/42 CEE EN 60601-1 and EN 60601-1-2.
0476	Notified body: CERMET.
	MET Mark
<small>UL 60601-1 CSA C22.2 No. 601-1 - E113015</small>	UL-CSA conformity
	Manufacturer
QTY.1	Do not allow fingers to contact moving parts
	Applied part of type "B" as per norm EN 60601-1
	Can be sterilized in autoclave up to a maximum temperature of 135° C
	Non-sterile
	The device and its accessories must not be disposed of or treated as solid urban wastes
I	Activation switch "on"
O	Activation switch "off"
	Alternating current
	Connection of the foot pedal
Rx Only	For US market only <b>CAUTION</b> US Federal Law restricts this device to sale by or on the order of a dentist

## 02 → IDENTIFICATION DATA


An exact description of the model and of the serial number of the device will enable our After-Sales Service to provide fast and efficient support.

Always refer these data whenever you contact an Authorized Mectron Service Center.

EN

### 02.1 → IDENTIFICATION PLATE OF THE DEVICE

Every device has an identification plate that bears the technical characteristics and the serial number. The identification plate is located on the inferior panel of the device. Additional specifications are reported in this manual (See Chapter 11 - TECHNICAL DATA).


### 02.2 → IDENTIFICATION DATA OF THE HANDPIECE

The handpiece is laser-marked with Mectron logo (ref. 1) and with its serial number (ref. 2).

### 02.3 → IDENTIFICATION DATA OF THE INSERTS

The following data are laser-marked on each insert: the name of the insert (ref. 3), the Mectron logo (ref. 4) and the lot number to which the insert belongs (ref. 5).

**03 → DELIVERY****03.1 → LIST OF THE COMPONENTS OF THE MULTIEZO**

See inside cover

EN

The multipiezo consists of:

- A** device body
- B** peristaltic pump
- C** scaler handpiece complete with LED light front cone
- D** metal front cone without light
- E** bottle complete with caps
- F** screw-on lids for the bottles
- G** inserts/insert kit
- H** mectron torque wrench
- I** electrical power supply cable
- L** foot pedal with cord and connector
- M** manual of use and maintenance

The multipiezo consists of accessories that can be ordered separately. Refer to the “Packing List” included in your package, so that you may exactly know the quantity and type of the accessories supplied with the device you have purchased.

The package of the device is sensitive to strong collisions, because it contains electronic components.

Therefore, special precautions must be taken for transport and storage.

Do not overlap multiple boxes, in order not to squash the packages underneath.

All the materials shipped by Mectron have been inspected upon their delivery.

The device is delivered duly protected and packed.

When receiving the device, check for the possible presence of damages incurred during the transport and, should that be the case, file a complaint with the transporter.

Preserve the package in the event that you need to send the device to an Authorized Mectron Service Center and so that you can place it back inside during prolonged periods of disuse.

 **WARNING:** Before starting to operate with the device, make sure that you have stock material (handpiece, inserts, wrenches) available to use in case the device stops due to a fault or of inconveniences.

## 04 → INSTALLATION

### 04.1 → FIRST INSTALLATION

The device must be installed in a comfortable place suitable for its use.

The multipiezo can be purchased ready for it to be used instantly or may have to be enabled by typing an activation key.

In the case that your device requires an activation key, the procedures to follow may vary from country to country.

Always refer to your retailer to obtain information on the above.

### 04.2 → SAFETY PRECAUTIONS DURING INSTALLATION

**⚠ WARNING: Contraindications. Interference with other equipment.** Though compliant with the standard IEC 60601-1-2, the multipiezo may nonetheless interfere with other devices nearby. The multipiezo must not be used near to or stacked on other devices. However, if this were to prove necessary, you must check and monitor correct operation of the device in that configuration.

**⚠ WARNING: Contraindications. Interference from other equipment.** An electrical scalpel or other electro-surgical units near the multipiezo device may interfere with its correct operation.

**⚠ CAUTION:** The electrical system of the premises in which the device is installed and used must be compliant to the norms in force and to the relative electrical safety precautions.

**⚠ CAUTION:** To avoid any risk of electric shock this device must be grounded.

**⚠ WARNING: Do not install the device in places where there is a risk of explosion.**

The device cannot operate in environments where anesthetic or flammable mixtures are present.

**⚠ WARNING:** Install the device in a place protected against collisions or against accidental sprays of water or liquids.

**⚠ WARNING:** Do not install the device above or near heat sources. Foresee adequate air circulation around the device when installing it.

**⚠ CAUTION:** Do not expose the device to direct sunlight or to sources of UV light.


**⚠ CAUTION:** The device can be transported, but it must be handled with care when it is displaced. Position the foot pedal on ground, so that it can only be activated intentionally by the operator.


**⚠ CAUTION:** Before connecting the handpiece to its cord, make sure that the electrical contacts are perfectly dry, on both sides. If need be, dry them with compressed air.

### 04.3 → CONNECTION OF THE ACCESSORIES

Connect the foot pedal to the back of the device in the socket marked with the symbol  by means of the plug of the pedal cable, until you hear a “click” sound.


1

2


Insert the power supply cable in its connection located at the back of the device. Connect it to the power outlet;

3

Unscrew the cap from the bottle and fill it in with the needed solution;


⚠ **CAUTION:** each bottle can contain liquids up to a maximum of 500 ml.


4

Check that the hose connected to the bottle's cap is correctly installed, then screw the cap back on its bottle;

⚠ **CAUTION:** check that the female coupling of the cap is clean and free from impediments.

⚠ **CAUTION:** check that the male coupling on the unit is clean and that its o-rings are not consumed.

5

Keep the bottle in a vertical position and push it towards the unit's body until it is firmly connected;

⚠ **CAUTION:** Do not kipsize the bottle as its cap is not watertight. The leaking of the physiological solution or of aggressive liquids can damage the surfaces.


## MULTIPIEZO

EN

Correctly connect the handpiece on its cord making sure to match the dot on the handpiece connector with the groove on the cord connector.


Make sure that the electrical contacts of both handpiece and cord are perfectly dry; if need be, dry them by blowing compressed air;

6

Position the handpiece on its support.

7

## 05 → USE

### 05.1 → SWITCHING THE DEVICE ON AND OFF

#### Turning the device on


With the device in view from its front side, position the switch on the left of the device body on position "I", being careful not to press the foot pedal.

4 symbols appear on the device (ref. U inside cover) that turn off one by one. At this point, the device sets on the default setting and is ready for it to be used.

#### Turning the device off

With the device in view from its front side, position the switch on the left of the device body on position "O", being careful not to press the foot pedal.

The device turns off.


**PLEASE NOTE** whenever the device is started, the following default setting is programmed:

- "power" 1
- "light" AUTO
- "function" ENDO
- "irrigation" 3

**⚠ CAUTION:** Always position the device in way so that the power switch is easily reachable, since this switch is considered as a load-break switch.

## 05.2 → DESCRIPTION OF THE KEYBOARD

### Touch keyboard

The user can configure the device by simply touching on the touch keyboard. Depending on the selected setting, the electronic feedback system automatically adjusts the correct operating frequency.

### POWER (ref. N inside cover)

The power can be adjusted by selecting the numbers on the touch keyboard in the “power” column. 7 power levels are foreseen: 0 = the handpiece has no power: only the irrigation outflows from the handpiece.

From 1 to 6= the power can be increased in an incremental way.

The power level can be regulated in each function, as indicated here below:

ENDO - 7 power levels: from 0 to 6

PERIO - 7 power levels: from 0 to 6

SCALER - 7 power levels: from 0 to 6

### FILLING AND FLUSHING THE IRRIGATION CIRCUIT (ref. O inside cover)

The device is equipped with the “PUMP/ FLUSH” key which, depending on the mode of use, allows to perform the PUMP function or the FLUSH function.

The PUMP function can be used at the beginning of the treatment, to flush the entire irrigation line up to the insert, so that the treatment can be started with the necessary irrigation (see Paragraph 05.4).

The FLUSH function (see Chapter 06 - FLUSH FUNCTION) allows to perform a cleaning cycle of the irrigation circuit.

This function must be performed when you have finished using the device and before cleaning, disinfecting, and sterilizing all the parts (see Chapter 08 - CLEANING AND STERILIZATION).

# MULTIPIEZO

EN



## LIGHT (ref. P inside cover)

Depending on the type of treatment that needs to be performed, it is possible to choose 3 options from the “light” list:

- By choosing the AUTO option, the LED light on the front terminal of the handpiece is switched on by pressing the foot pedal, and automatically switches off 3 seconds after the pedal is released;
- By choosing the ON option, the LED light on the front terminal of the handpiece stays on permanently, regardless of whether pressure is applied on the foot pedal or not. The light switches off 100 seconds after the last pressure of the foot pedal and the option shifts from ON to AUTO;
- By choosing the OFF option, the LED light on the front terminal of the handpiece stays off permanently.

The position of the LED light on the front terminal of the handpiece can be adjusted in the following way:

- Hold the body of the handpiece and lightly unscrew the metal ring nut located at the base of the front terminal, rotating it counter-clockwise;
- Rotate the front terminal so that the LED light goes into the desired and necessary position;
- To fasten it into position, screw the metal ring nut, rotating it clockwise.


## FUNCTIONS (ref. Q inside cover)

Depending on the type of treatment it is possible to choose one of the 3 options available from the “function” list, as follows:

**ENDO:** dedicated to endodontic treatments as cleaning of root canals and retro approach.

**PERIO:** dedicated to the supra- and sub-gingival prophylaxis procedures, to root planing and to the cleaning of the implant surface.

**SCALER:** dedicated to scaling, restoration and prosthesis techniques.


### SOFT MODE (ref. R inside cover)

Once the PERIO or the SCALER function is selected, you can make the treatment more delicate by choosing the “soft mode”.

The function is not available in ENDO.  
Indicated for very sensitive patients.

### IRRIGATION (ref. S inside cover)

The capacity of the peristaltic pump can be adjusted by selecting the numbers on the touch keyboard in the “irrigation” column. 7 capacity levels are foreseen:

0 = operation of the pump is closed: no irrigation outflows from the insert.


From 1 to 6= the pump flow goes from 10 ml/min to approximately 40 ml/min.

The possibility of choosing the irrigation capacity levels is tied to the type of function selected, in the following manner:

ENDO - 7 flow levels: from 0 to 6

PERIO - 7 flow levels: from 0 to 6.

SCALER - 7 flow levels: from 0 to 6.

### SYMBOLS (ref. T inside cover)

The multipiezo is equipped with a diagnostic circuit that allows to detect operating abnormalities and to view their type on the keyboard via their relative symbol.

To help the user identify the malfunctioning part, four symbols are foreseen which are described in paragraph 12.1.

## 05.3 → SAFETY PRECAUTIONS BEFORE AND DURING USE

EN

**⚠ WARNING:** Before starting to operate with the device, make sure that you have stock material (handpiece, inserts, wrenches) available to use in case the device stops due to a fault or of inconveniences.

**⚠ WARNING:** Only use original Mectron inserts, accessories, and spare parts.

**⚠ WARNING: Use of non-original Mectron inserts:** this use entails finite damage to the handpiece threading, thus compromising correct operation and risking to cause harm to the patient.

**⚠ CAUTION: Contraindication. Ultrasonic scaler.** Do not perform treatments on prosthetic artifacts made of metal or ceramics. The ultrasonic vibrations could lead to the de-cementing of the artifacts.

**⚠ WARNING: Contraindications.** Do not use the multipiezo on patients who carry heart stimulators (Pace-makers) or other implantable electronic devices. This precaution also applies to the operator.

**⚠ WARNING: Contraindications.** Do not perform scaling treatments without water spray in order to avoid the insert overheating which may cause damages to the tooth. Treatments with no water spray can be carried out only with "Dry Work" inserts which do not have the water passage.

**⚠ CAUTION:** Only use inserts with water passage when performing treatments that require irrigation.

**⚠ WARNING: Treatments that require irrigation.** Always check operation of the irrigation before and during use. Make sure the fluid outflows from the insert. Do not use the device if the irrigation does not work or if the pump is defective.

**⚠ WARNING: Checking device status before the treatment.** Always check that there is no water underneath the device. Before every treatment, always check that the device works perfectly and that

the accessories are efficient. In case you uncover operating abnormalities, do not perform the treatment. Contact an Authorized Mectron Service Center if the abnormalities concern the device.

**⚠ WARNING: Infections control.**

**First use:** All new and repaired accessories are supplied in NON STERILE conditions. Before use, and after each treatment, they must be cleaned and sterilised in strict compliance with the instructions given in Chapter 08 - CLEANING AND STERILIZATION.

**Subsequent uses:** After every treatment, clean and sterilize all the reusable parts and accessories, following the instructions provided in Chapter 08 - CLEANING AND STERILIZATION.

**⚠ WARNING:** To ensure that the handpiece cools down, always activate it with the irrigation circuit correctly installed and filled. To fill the irrigation circuit, always use the FLUSH function.

**⚠ CAUTION:** To correctly use the device, you must press the foot pedal and start it with the insert not making contact with the part to be treated, so that the electronic circuit can recognize the best point of resonance of the insert without interferences, allowing its optimal performance.

**⚠ WARNING:** Before every treatment, make sure that the insert appropriate for the treatment is inserted on the handpiece. Exclusively use the Mectron torque wrench to fasten the insert to the handpiece.

**⚠ WARNING:** The patient must not come into contact with the device body or the foot pedal.

**⚠ WARNING:** Do not change the insert while the handpiece is operating, to prevent causing wounds to the operator.

**⚠ CAUTION: FLUSH function.** After the device is used with aggressive and non-aggressive solutions, it is necessary to

perform a cleaning cycle on the tubes and the handpiece with the FLUSH function (see Chapter 06 - FLUSH FUNCTION). If the tubes are not cleaned, the crystallization of the salts may seriously damage the device.

⚠ **CAUTION: FLUSH function.** The FLUSH function must be used after every treatment, before starting the cleaning and sterilization procedures.

⚠ **WARNING: Breakage and wear-out of the inserts.** High frequency oscillations and wear-out may, in rare circumstances, lead to the breakage of the insert.

Do not bend, change shape of, or re-sharpen an insert in any way.

Bending an insert or applying leverage on it can lead to its breakage.

Deformed or otherwise damaged inserts are susceptible to breakage during their use. These inserts must never be used.

Excessive pressure on the inserts during their use can lead to their breakage.

If an insert breaks, check that none of its fragments remain in the treated part and, at the same time, apply effective suction to remove them.

The patient must be instructed to breathe through his nose during the treatment, or a dental dam must be used to prevent

the patient from ingesting fragments of broken inserts.

When the nitride coating wears out, the cutting efficiency decreases; re-sharpening the insert damages it and is therefore forbidden. Check that the insert is not worn out.

During the intervention, frequently check that the insert is intact, especially its top.

During the intervention, avoid prolonged contact with retractors or with metallic instrumentation in use.

⚠ **CAUTION: Contraindications.** Allow reusable, autoclavable items (the handpiece, the inserts, the torque wrench, and any other accessory that can be sterilized ) to gradually return to room temperature after steam sterilization and prior to usage. The cooling process must not be accelerated.

⚠ **CAUTION: The electrical contacts of both handpiece and cord connectors must be dry.** Before connecting the handpiece to its cord make sure that the electrical contacts of both connectors, are perfectly dry, especially after the sterilization cycle in autoclave. If need be, dry the contacts by blowing compressed air onto them.

## 05.4 → INSTRUCTIONS FOR USE

After having connected all the accessories as illustrated in Paragraph 04.3 proceed as follows:

EN


To fill the irrigation circuit use the PUMP function by selecting PUMP/FLUSH on the touch keyboard: all the other selection options present on the display are disabled and the PUMP/FLUSH writing flashes;

1


While the PUMP/FLUSH writing flashes, press the foot pedal once and release it, the PUMP/FLUSH stops flashing and the irrigation circuit begins to be filled up;

2


As soon as the peristaltic pump starts, the entire scale of values of the "irrigation" section lights up and during the liquid passage, the value of the irrigation shifts from 6 to 0;


3

The cycle can be interrupted as soon you see liquid outflowing from the multipiezo handpiece, by pressing on PUMP/FLUSH menu or, alternatively, by pressing the foot pedal. The function PUMP is disabled and the keyboard is enabled again, and displays the last setting used;


4

5


Screw the chosen insert onto the multi-piezo handpiece, until it reaches the limit stop;

6


Tighten the insert by using the Mectron torque wrench.

To correctly use the Mectron torque wrench, operate as follows:

7

Put the insert inside the wrench as shown;

8

Firmly hold the central body of the handpiece;

**!** **CAUTION:** The handpiece must not be grabbed by its terminal part and/or cord, but only by its central body.

The handpiece must not be rotated, but must be grasped firmly, and you must only rotate the wrench.


Rotate the wrench clockwise until the friction snaps (the external body of the wrench rotates compared to the handpiece, emitting mechanical "CLICK" sounds). The insert is now perfectly tightened;

## MULTIPIEZO

On the keyboard, select the type of function and irrigation necessary and the light, if desired;


Lift the handpiece and press the foot pedal in order to start the treatment;

9

At the end of the treatment position the scaler handpiece on its support.

10

## 05.5 → IMPORTANT INFORMATION ON THE INSERTS

### ⚠ WARNING:

- When the titanium nitride coating is visibly worn out, the insert must be replaced. Use of an overly worn out insert reduces its efficiency.
- Diamond-coated inserts: the diamond inserts must be replaced when the titanium nitride coating is visibly worn out and, at any rate, after a maximum of 10 treatments.
- Do not activate the handpiece while the insert is in contact with the treated part, so that the electronic circuit can recognize the best point of resonance of the insert and allow its optimal performance.
- Check the condition of wear of the insert and that it is intact before and during every use. In the event that a drop in performance becomes apparent, proceed to replace it.
- Use original Mectron inserts only. Use of non-original inserts, in addition to voiding the warranty, damages the threading of the multipiezo handpiece, with the risk of no longer being able to screw the original inserts correctly during subsequent use. Moreover, the device settings are tested and guaranteed to operate correctly only when original Mectron inserts are used.
- Do not change the shape of the insert in any way by either bending or filing it. This could cause it to break.
- Do not use an insert that has suffered any type of deformation.
- Do not attempt to sharpen the insert used.
- Always check that the threaded parts of the insert and of the handpiece are perfectly clean – see Chapter 08 - CLEANING AND STERILIZATION.
- If excessive pressure is applied on the insert, it can cause the insert to break and possibly harm the patient.


## 06 → FLUSH FUNCTION

The FLUSH function allows to run a cleaning cycle on the irrigation circuit of the handpiece used during the treatment, by following the steps described in this chapter.


⚠ **CAUTION: FLUSH function.** The FLUSH function must be used after every treatment, before starting the cleaning and sterilization procedures.

⚠ **CAUTION:** If the tubes are not cleaned, the crystallization of the salts may seriously damage the device.

EN

- 1 Fill in the bottle with water, possibly demineralized water.  
Connect the bottle to the device;


- 2 Position the handpiece, with or without insert, above a container to contain the liquid that will outflow during the cleaning cycle;

- 3 To enter the flushing mode, select PUMP/FLUSH menu on the touch keyboard: all the other selection options present on the display are disabled;

**PLEASE NOTE** You can exit the flushing mode at any time by pressing PUMP/FLUSH menu, the keyboard is enabled again, and displays the last setting used;

- 4 To make the FLUSH function start while the PUMP/FLUSH writing is flashing, press the foot pedal once and release it: the screen stops flashing and the cleaning cycle starts;


**PLEASE NOTE:** the FLUSH function can be interrupted at any time by pressing the PUMP/FLUSH writing again or, alternatively, by pressing on foot pedal. The keyboard is enabled again, and displays the last setting used;

## MULTIPIEZO

EN

As soon as the peristaltic pump starts, the entire scale of values of the “irrigation” section lights up and during the liquid passage, the value of the irrigation shifts from 6 a 0 . The cycle lasts 20 seconds. Once it has ended, the keyboard is enabled again, and displays the last setting used.

5

After the cleaning cycle of the irrigation circuit move on to disassembling the different parts (see Chapter 07 - DISASSEMBLY OF PARTS FOR


CLEANING AND STERILIZATION) and to their cleaning and sterilization (see Chapter 08 - CLEANING AND STERILIZATION).

## 07 → DISASSEMBLY OF PARTS FOR CLEANING AND STERILIZATION


Before moving on to the cleaning procedures described in Chapter 08 - CLEANING AND STERILIZATION, disconnect all the accessories and components of the multipiezo.

**⚠ WARNING: Switch the device off.**  
Always turn the device off via the switch, and disconnect the power supply cable from the power outlet and from the device body, before performing the cleaning and sterilization interventions.


1

2


If present, unscrew the insert from the handpiece by using the torque wrench;


3

Disconnect the handpiece from its cord;


⚠ **CAUTION:** Do not attempt to unscrew or turn the connector when disconnecting the handpiece. The connector could get damaged.

4


Unscrew the front cone of the handpiece;

NOTE: on the cone with LED light the metal ring cannot be separated from the cone in plastic;

5

Disconnect the bottle from the device body by pulling them outwards;

6

Disconnect the foot pedal from the device: Grab the connector of the foot pedal, press the release flap, and pull the connector back;

⚠ **CAUTION:** Do not attempt to unscrew or turn the connector during the disconnection: the connector could get damaged.

⚠ **CAUTION:** When disconnecting the foot pedal, always and only hold the connector of the cord. Never actually pull on the cord itself.

## 08 → CLEANING AND STERILIZATION

### 08.1 → CLEANING OF THE FOOT PEDAL

EN

#### PREPARATION

- Disconnect the foot pedal from the device body (see Chapter 07)

- ⚠ **CAUTION:** Do not sterilize the foot pedal, it would stop working.
- ⚠ **CAUTION:** Do not clean the foot pedal under running tap water.
- ⚠ **CAUTION:** Do not immerse the foot pedal in liquids and/or solutions of various nature.

#### MATERIALS NECESSARY

- Clean, soft cloths with low fiber-release
- Detergent solution (pH 6-9) and, if need be, non-aggressive disinfectant solution with neutral pH (pH7)


- ⚠ **CAUTION:** We recommend that you use water-based disinfectant solutions with neutral pH (pH7). Alcohol-based disinfectant solutions and hydrogen peroxide are contra-indicated, because they can fade the color and/or damage the plastic materials. This also holds true for chemical products such as acetone and alcohol.

#### CLEANING AND DISINFECTION METHOD - Foot pedal


Clean the surface of the foot pedal with a clean, soft cloth with low fiber-release, dampened with a detergent solution (pH 6-9) and, if need be, disinfect with a non-aggressive disinfectant solution with neutral pH (pH7), following the instructions provided by the manufacturer of the solution;

Dry the foot pedal with a clean, non-abrasive cloth with low fiber-release .

1

2

## 08.2 → CLEANING OF THE DEVICE BODY

## PREPARATION

- Make sure that all the following accessories have been removed/disconnected from the device body (see Chapter 07):
  - electrical power supply cable
  - foot pedal
  - handpiece
  - bottle

**⚠ WARNING:** Do not sterilize the device body: it could stop working and cause harm to person and/or damages to things.

**⚠ WARNING:** Always turn the device off via the switch and disconnect it from the electrical network before performing the cleaning and disinfection interventions.


**⚠ WARNING:** The device is not protected against penetration of liquids. Do not spray liquids directly onto the surface of the device.

## MATERIALS NECESSARY

- Clean, soft cloths with low fiber-release
- Detergent solution (pH 6-9) and, if need be, non-aggressive disinfectant solution with neutral pH (pH7)


**ⓘ CAUTION:** We recommend that you use water-based disinfectant solutions with neutral pH (pH7). Alcohol-based disinfectant solutions and hydrogen peroxide are contraindicated, because they can fade the color and/or damage the plastic materials. This also holds true for chemical products such as acetone and alcohol.

## CLEANING AND DISINFECTION METHOD - Device body

1

Clean the surface of the device and of the handpiece cord with a clean, soft cloth with low fiber-release, dampened with a detergent solution (pH 6-9) and, if need be, disinfect with a non-aggressive disinfectant solution with neutral pH (pH7), following the instructions provided by the manufacturer of the solution;

2

Dry the device body and the handpiece cord with a clean, non-abrasive cloth with low fiber-release.

## 08.3 → CLEANING AND STERILIZATION OF THE HANDPIECE

### PREPARATION

- Perform the FLUSH function (see Chapter 06)
- Disconnect the handpiece from the cord attached to the device body (see Chapter 07)
- If present, disconnect the insert from the handpiece with the torque wrench (see Chapter 07)
- Unscrew the front cone (see Chapter 07)

⚠ **CAUTION:** Do not attempt to unscrew or turn the connector when disconnecting the handpiece. The connector could get damaged.

⚠ **CAUTION:** Do not immerse the handpiece in disinfectant solutions or other liquids, because it could get damaged.

⚠ **CAUTION:** Do not immerse the handpiece in an ultrasonic tank.

### MATERIALS NECESSARY


- Enzymatic detergent with pH 6-9
- Water
- Container for solution with enzymatic liquid
- Clean and soft cloth with low fiber-release
- If need be, non-aggressive disinfectant solution with neutral pH (pH7)
- Brush with soft nylon bristles
- Distilled water
- Compressed air
- Single-use sterilization packets
- Steam sterilization


⚠ **CAUTION:** If you intend to disinfect, we recommend that you use water-based disinfectant solutions with a neutral pH (pH7). Alcohol-based disinfectant solutions and hydrogen peroxide are contraindicated, because they can fade the color and/or damage the plastic materials. This also holds true for chemical products such as acetone and alcohol. Always rinse with sterile water to preserve the disinfection.

⚠ **CAUTION:** Once used, correctly dispose of the enzymatic detergent solution, do not recycle.

### CLEANING AND DISINFECTION METHOD - Handpiece


Prepare an enzymatic detergent solution with pH 6-9, as per the instructions of the manufacturer;

2



Clean the surface of the handpiece and of its connector with a clean, soft cloth with low fiber-release, dampened with a detergent solution (pH 6-9) and, if need be, disinfect with a non-aggressive disinfectant solution with neutral pH (pH7), following the instructions provided by the manufacturer of the solution;

3


Delicately clean the surface of the handpiece using the enzymatic detergent solution with the aid of a brush with soft nylon bristles, placing special care to the following areas:

- threading of the handpiece
- titanium stem
- front terminal in its external and internal parts

4

Accurately rinse with tap water to eliminate any detergent residues and perform the last rinse with distilled water;

5

Thoroughly dry all the parts, especially the electrical contacts, blowing on them with compressed air;


## MULTIPIEZO

EN


Once you've finished the cleaning operations, check all the objects under an adequate source of light, paying attention to the details that could conceal dirt residues (threading, cavities, grooves) and, if need be, repeat the cleaning cycle. Finally, check the integrity of those parts and elements that may have worn out during their use;

Individually seal the handpiece (without inserts) and the front terminal, separately, in single-use packets for sterilization. Proceed with the sterilization.

6

7

### STERILIZATION METHOD - Handpiece

The handpiece is built with materials that resist to a maximum temperature of 135°C for a maximum time of 20 minutes.

The sterilization process in steam autoclave must be performed setting the parameters indicated below, which are validated by Mectron and reach the sterilization safety levels (SAL 10<sup>-6</sup>):

- 3 times Pre-vacuum (minimum pressure 60 mBar).
- Sterilization temperature 132°C (interval 0°C ÷ +3°C).
- Sterilization time 4 minutes.
- Minimum drying time 20 minutes.

All the sterilization phases must be performed by the operator in compliance to the norms UNI EN ISO 17665-1:2007, UNI EN ISO 556-1:2002 e ANSI/AAMI ST46:2002.

⚠ **CAUTION:** Do not sterilize the handpiece with the insert screwed onto it.


⚠ **WARNING: Checking infections - Parts that can be sterilized** - Diligently remove all residues of organic dirt before the sterilization.

⚠ **CAUTION:** Perform the sterilization using a water steam autoclave only. Do not use any other sterilization procedure (dry heat, irradiation, ethylene oxide, gas, low temperature plasma, etc.).

⚠ **CAUTION:** Do not exceed the allowed load of the steam sterilizer.

⚠ **CAUTION:** At the end of the sterilization cycle, let the handpiece cool down completely prior to using it.

⚠ **CAUTION:** The electrical contacts of the connector must be dry. At the end of the sterilization cycle and before connecting the handpiece to the cord, make sure that the electrical contacts of their connectors are perfectly dry; if need be, dry them by blowing compressed air.

## 08.4 → CLEANING AND STERILIZATION OF THE INSERTS

## PREPARATION

- Disconnect the insert from the handpiece with the torque wrench (see Chapter 03)

⚠ **CAUTION:** Always disconnect the insert from the handpiece before proceeding to its cleaning and sterilization.

EN


## MATERIALS NECESSARY

- Enzymatic detergent with pH 6-9
- Water
- Ultrasound tank
- Brush with soft nylon bristles
- Compressed air
- Distilled water
- Single-use sterilization packets
- Steam sterilization
- Syringe

⚠ **CAUTION:** Do not use hydrogen peroxide. If you intend to disinfect the inserts, do not use hydrogen peroxide, but only disinfectants with neutral pH (pH7); always rinse with sterile water to preserve the disinfection.


⚠ **CAUTION:** Once used, correctly dispose of the enzymatic detergent solution, do not recycle.

## CLEANING AND DISINFECTION METHOD - Inserts

1

Prepare an enzymatic detergent solution with pH 6-9, as per the instructions of the manufacturer;

2


Place the insert in the ultrasound tank submerged by the enzymatic detergent solution at 40°C, for at least 10 minutes;

## MULTIPIEZO

EN


Remove the insert from the ultrasonic tank and rinse with distilled water;

3


Delicately brush the surface of the insert again with the brush with soft nylon bristles;

4


Rinse the internal canal of the insert with distilled water injected under pressure (3,8 bar) for at least 10 seconds to eliminate all residues;


5

Once you've finished the cleaning operations, check all the inserts under an adequate source of light, paying attention to the details that could conceal dirt residues (threading, cavities, grooves) and, if need be, repeat the cleaning cycle;


6

7

⚠ **CAUTION:** Before starting the sterilization cycle, make sure that the insert is thoroughly dry, both externally and internally. For this purpose, blow compressed air both externally and through the internal thru hole; this will prevent the onset of stains and haloes on the surface or rust inside the insert.

8

Individually seal the inserts in single-use sterilization packets. Sterilize the inserts in a steam autoclave.

EN

### STERILIZATION METHOD - Inserts

The sterilization process in steam autoclave must be performed setting the parameters indicated below, which are validated by Mectron and reach the sterilization safety levels (SAL 10<sup>-6</sup>):

- 3 times Pre-vacuum (minimum pressure 60 mBar).
- Sterilization temperature 132°C (interval 0°C ÷ +3°C).
- Sterilization time 4 minutes.
- Minimum drying time 20 minutes.

All the sterilization phases must be performed by the operator in compliance to the norms UNI EN ISO 17665-1:2007, UNI EN ISO 556-1:2002 e ANSI/AAMI ST 46:2002.

⚠ **CAUTION:** Do not sterilize the insert screwed onto the handpiece.

⚠ **WARNING: Checking infections - Parts that can be sterilized** - Diligently remove all residues of organic dirt before the sterilization.

⚠ **CAUTION:** Perform the sterilization using a water steam autoclave only. Do not use any other sterilization procedure (dry heat, irradiation, ethylene oxide, gas, low temperature plasma, etc.).

⚠ **CAUTION:** Do not exceed the allowed load of the steam sterilizer.

## 08.5 → CLEANING AND STERILIZATION OF THE TORQUE WRENCH

### PREPARATION

Take the wrench

EN

### MATERIALS NECESSARY

- Enzymatic detergent with pH 6-9
- Water
- Ultrasound tank
- Brush with soft nylon bristles
- Distilled water
- Clean and soft cloth with low fiber-release
- Medical degree lubricant
- Single-use sterilization packets
- Steam sterilization


⚠ **CAUTION: Do not use hydrogen peroxide.** If you intend to disinfect the wrench, do not use hydrogen peroxide, but only disinfectants with neutral pH (pH7); always rinse with sterile water to preserve the disinfection.

⚠ **CAUTION: Once used, correctly dispose of the enzymatic detergent solution, do not recycle.**

### CLEANING AND DISINFECTION METHOD - Torque wrench


Prepare an enzymatic detergent solution with pH 6-9, as per the instructions of the manufacturer;


1

Place the insert in the ultrasound tank submerged by the enzymatic detergent solution at 40°C, for at least 10 minutes;


2


3

Rinse with tap water to eliminate all residues;


4

Delicately brush the surface of the wrench with a brush with soft nylon bristles using distilled water;

5


Dry the wrench with a soft cloth with low fiber-release;

6

Lubricate with medical degree lubricant in the indicated point;

**⚠ CAUTION:** do not use lubricants based on oil or silicon.


7

Once you've finished the cleaning operations, check the wrench under an adequate source of light, paying attention to the details that could conceal dirt residues (cavities, grooves) and, if need be, repeat the cleaning cycle;

## MULTIPIEZO

Seal the wrench individually in a single-use sterilization packet.  
Sterilize the wrench in a steam autoclave.

8

EN

### STERILIZATION METHOD - Torque wrench

The sterilization process in steam autoclave must be performed setting the parameters indicated below, which are validated by Mectron and reach the sterilization safety levels (SAL 10<sup>-6</sup>):

- 3 times Pre-vacuum (minimum pressure 60 mBar).
- Sterilization temperature 132°C (interval 0°C ÷ +3°C).
- Sterilization time 4 minutes.
- Minimum drying time 10 minutes.

All the sterilization phases must be performed by the operator in compliance to the norms UNI EN ISO 17665-1:2007, UNI EN ISO 556-1:2002 e ANSI/AAMI ST 46:2002.

**⚠ WARNING: Checking infections - Parts that can be sterilized** - Diligently remove all residues of organic dirt before the sterilization.

**⚠ CAUTION:** Perform the sterilization using a water steam autoclave only. Do not use any other sterilization procedure (dry heat, irradiation, ethylene oxide, gas, low temperature plasma, etc.).

**⚠ CAUTION:** do not exceed the allowed load of the steam sterilizer.

## 08.6 🔄 CLEANING OF THE BOTTLE

### PREPARATION

- Disconnect the bottle from the device body (see Chapter 07);
- Unscrew the cap from the bottle.


**⚠ CAUTION:** Do not sterilize the bottle and its cap in autoclave, they could get damaged.

### MATERIALS NECESSARY

- Water
- Detergent solution (pH 6-9). If need be, non-aggressive disinfectant solution with neutral pH (pH7)
- Clean and soft cloth with low fiber-release
- Sterile water

**⚠ CAUTION:** If you intend to disinfect the bottle and its cap we recommend that you use water-based disinfectant solutions with a neutral pH (pH7). Alcohol-based disinfectant solutions and hydrogen peroxide are contra-indicated, because they can fade the color and/or damage the plastic materials. This also holds true for chemical products such as acetone and alcohol. Always rinse with sterile water to preserve the disinfection.

## CLEANING AND DISINFECTION METHOD - Bottle


1

Carefully rinse under running water the internal and the external surfaces of the bottle and of the cap;


2

Clean the external surfaces of the bottle and of the cap with a clean, soft cloth with low fiber-release, dampened with a detergent solution (pH 6-9) and, if need be, disinfect with a non-aggressive disinfectant solution with neutral pH (pH7), following the instructions provided by the manufacturer of the solution;

3

Remove any detergent remains with a soft cloth with low fiber-release dampened with sterile water in order to maintain the disinfection;

4

Once you've finished the cleaning operations, check the bottle and the cap under an adequate source of light, paying attention to the details that could conceal dirt residues (cavities, grooves) and, if need be, repeat the cleaning cycle.

## 09 → MAINTENANCE

EN

If the device is not used for prolonged time, observe the following recommendations:

- 1 Run a complete cleaning cycle on the irrigation circuit with the FLUSH function (see Chapter 06 - FLUSH FUNCTION);
- 2 Empty the water circuits by disconnecting the bottle and by letting the handpiece function for some seconds;
- 3 Disconnect the device from the electrical network;
- 4 If the period of disuse is prolonged, put the device back in its original package and store it in a safe place;

- 5 Prior to using the device again, clean and sterilize the handpiece, the inserts and the wrench following the instructions provided in the Chapter 08 - CLEANING AND STERILIZATION;
- 6 Check that the inserts are not worn out, deformed, or broken, placing special attention to the integrity of their tip.

**⚠ WARNING:** Periodically check that the electrical power cable is intact; if it is damaged, replace it with an original Mectron spare part.


### 09.1 → REPLACEMENT OF THE PERISTALTIC PUMP

The peristaltic pump on the back of the device is protected by a plastic cover. This cover can be removed by pressing its sides and pulling them backwards;


**⚠ CAUTION:** Before doing any kind of operation on the peristaltic pump, make sure that the device is disconnected from the power outlet and that the liquid container is not connected.


Disconnect the two silicon tubings of the peristaltic pump from their respective couplings, which are positioned beneath the pump;

1


2


3

Take out the peristaltic pump from its seat by pulling it backwards;

4

Connect the new peristaltic pump to its seat until a “click” sound is heard and then connect the two pump tubings to their respective couplings positioned beneath the pump;

5

Re-position the plastic cover on the peristaltic pump.

## 10 → DISPOSAL MODES AND PRECAUTIONS

### ⚠ **WARNING: Hospital wastes.**

Hospital wastes. Treat the following items as hospital waste:

- Inserts: when they are worn out or broken;
- Torque wrench for inserts: when worn out or broken.

The use and throw-away materials and materials that entail a biological risk must be disposed of in accordance to the norms in force locally concerning hospital wastes.

The multipiezo must be disposed of and treated as a waste for separate collection. Disregard of the previous points may entail a fine, pursuant to Directive 2002/96/CE. It is up to the purchaser to hand over the device for its disposal to the retailer who supplies him new equipment; the instructions for proper disposal are available from Mectron.

## 11 → TECHNICAL DATA

EN

<b>Device compliant to Dir. 93/42/CEE:</b>	Class IIa
<b>Classification as per EN 60601-1:</b>	I Applied part type B (handpiece, insert) IP 20 (device) IP 22 (foot pedal)
<b>Device for intermittent operation:</b>	55sec. ON - 30sec. OFF with irrigation 30sec. ON - 120sec. OFF without irrigation (ENDO, PERIO, SCALER)
<b>Power supply voltage:</b>	100-240 V ~ 50/60 Hz
<b>Max. power absorbed:</b>	90 VA
<b>Fuses:</b>	Type 5 x 20 mm T 2AL, 250V
<b>Operating frequency:</b>	Automatic scan From 24 KHz to 36 KHz
<b>Power types:</b>	ENDO PERIO SCALER SOFT MODE
<b>Peristaltic pump capacity:</b>	Adjustable on the touch screen 7 flow levels: from 0 (0 ml/min) to 6 (circa 28 ml/min) (see Paragraph 05.2 - IRRIGATION)
<b>LED system of the handpiece:</b>	Light function on AUTO: LED system of the handpiece switches on as soon as the device starts to operate and switches off 3 seconds after the foot pedal is released  Light function on ON: The LED of the handpiece is always on; after 100 seconds of foot pedal disuse it switches off by itself and the light function moves into position AUTO  Light Function on OFF: The LED of the handpiece is always off.
<b>Protections of the APC circuit:</b>	No handpiece detected Cord interruption Insert not tightened correctly or broken
<b>Operating conditions:</b>	from +10°C to +35°C Relative humidity from 30% to 75% Pressure of air P: 800hPa/1060hPa
<b>Transport and storage conditions:</b>	from -10°C to +70°C Relative humidity from 10% to 90% Pressure of air P: 500hPa/1060hPa
<b>Weights and sizes:</b>	2,4 Kg L - W - H 320 x 230 x 145 mm

## 11.1 → ELECTROMAGNETIC COMPATIBILITY EN 60601-1-2

**⚠ WARNING: Interference with other equipment.** Though compliant with the standard IEC 60601-1-2, the multipiezo may nonetheless interfere with other devices nearby.

The multipiezo must not be used near to or stacked on other devices. However, if this were to prove necessary, you must check and monitor correct operation of the device in that configuration.

**⚠ WARNING:** Portable and mobile radio communication appliances may affect the correct functioning of the device.

**⚠ WARNING: Interference from other equipment.** An electrical scalpel or other electro-surgical units near the multipiezo device may interfere with its correct operation.

**⚠ WARNING:** The device requires specific EMC precautions and must be installed and activated in accordance with the EMC information given in this paragraph.

## Guidance and manufacturer's declaration - Electromagnetic emissions

The multipiezo is intended for use in the electromagnetic environment specified below. The customer or user of the multipiezo should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – Guidance
RF emissions CISPR 11	Group 1	The multipiezo only uses RF energy for internal function. Therefore, its FR emissions are very low and are not likely to cause any interference with nearby electronic equipment.
RF emissions CISPR 11	Class B	The multipiezo is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

## Guidance and manufacturer's declaration - Electromagnetic immunity

The multipiezo is intended for use in the electromagnetic environment specified below. The customer or user of the multipiezo should ensure that it is used in such an environment.

Immunity test	Test level IEC 60601	Compliance level	Electromagnetic environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	The device continues to work regularly and in safety	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	The device continues to work regularly and in safety	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	The device continues to work regularly and in safety	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0,5 cycles 40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles 70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 s	The device can vary from the required levels of immunity with a duration of <5% / >95% / 5 s as long as the device remains in safety, no malfunctions have been detected and can be restored to pre-test status with the intervention of the operator	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	The device continues to work regularly and in safety	Power frequency magnetic fields should be at levels characteristic of a typical location in a commercial or hospital environment.

N.B.:  $U_T$  is the AC mains voltage prior to application of the test level.

**Guidance and manufacturer's declaration - Electromagnetic immunity**

The multipiezo is intended for use in the electromagnetic environment specified below. The customer or user of the multipiezo should ensure that it is used in such an environment.

EN

Immunity test	Test level IEC 60601	Compliance level	Electromagnetic environment Guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Veff from 150 kHz to 80 MHz</p> <p>3 V/m from 80 MHz to 2,5 GHz</p>	<p>The device continues to work regularly and in safety</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the device including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance <math>d = 1,2 \sqrt{P}</math></p> <p><math>d = 1,2 \sqrt{P}</math> from 80 MHz to 800 MHz  <math>d = 2,3 \sqrt{P}</math> from 800 MHz to 2,5 GHz</p> <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup>, may be less than the compliance level in each frequency range<sup>b</sup>.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

**N.B.:**

- (1) at 80 MHz and 800 MHz, the higher frequency range applies.
- (2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
- a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with any reasonable accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the multipiezo is used exceeds the applicable RF compliance level given above, the multipiezo should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the multipiezo.
- b Over the frequency range from 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

## Recommended separation distances between portable and mobile RF communications equipment and the multipiezo

The multipiezo is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the multipiezo can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile FR communications equipment (transmitters) and the multipiezo as recommended below, according to the maximum output power of the communications equipment.

EN

Rated maximum output power of transmitter 'W'	Separation distance according to the frequency of transmitter 'm'		
	from 150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	from 80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	from 800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be calculated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**N.B.:**

- (1) at 80 MHz and 800 MHz, the higher frequency range applies.
- (2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## 12 → TROUBLESHOOTING

### 12.1 → DIAGNOSTIC SYSTEM AND SYMBOLS ON KEYBOARD

The multipiezo is equipped with a diagnostic circuit that allows to detect operating abnormalities and to view their type on the keyboard via their relative symbol. By using the following chart, the user is guided toward the identification and possible solution of the malfunction detected.

EN

Symbols on keyboard	Possible cause	Solution
	Handpiece/cord contacts wet	Thoroughly dry the contacts with compressed air
	Handpiece not connected to device	Connect the handpiece
	Handpiece failure	Replace the handpiece
	Sync circuit malfunction	Contact an Authorized Mectron Service Center
	Insert not tightened correctly on handpiece	Unscrew the insert and correctly screw it on again with the Mectron torque wrench (see paragraph 05.4)
	Insert broken, worn-out or deformed	Replace the insert
	Handpiece/cord contacts wet	Thoroughly dry the contacts with compressed air
	Peristaltic pump malfunction	Check that there are no impediments to pump rotation Check that the pump and its two tubes are correctly installed
	The device has been turned off and on again without waiting 5 seconds	Turn device off and wait 5 seconds before turning it on again
	Abnormalities on electrical network or excessive electrostatic discharges or internal abnormalities	Turn device off and wait 5 seconds before turning it on again If the signal persists, contact an Authorized Mectron Service Center
	Turn-on procedure incorrect: the device has been turned on with the foot pedal pressed	Check that the foot pedal is not pressed. If the problem persists, disconnect the pedal and, if need be, contact an Authorized Mectron Service Center

12.2 → QUICK SOLUTION TO PROBLEMS

EN

Problem	Possible cause	Solution
<b>The device does not turn on after having brought the switch into position “I”</b>	The electrical power cable terminal is badly inserted in the rear plug of the device	Check that the power supply cable is firmly connected
	The electrical power cable is defective	Check that the power supply socket works properly Replace the electrical power cable
	The fuses are out of order	Replace the fuses (see paragraph 12.3)
<b>The device is on but not working. The display does not signal any error</b>	The foot pedal plug is incorrectly inserted in the device socket	Correctly insert the pedal plug in the socket on the back of the device
	The foot pedal does not work	Contact an Authorized Mectron Service Center
<b>The device is on but not working. One of the following symbols appears on the screen:</b> 	See paragraph 12.1 for the possible cause, according to the symbol that has been displayed	See paragraph 12.1 for the action to undertake, according to the symbol that has been displayed
<b>A slight whistling sound coming from the multipiezo handpiece is heard during operation.</b>	The insert is not correctly tightened on the handpiece	Unscrew and correctly screw the insert again with the Mectron torque wrench (see paragraph 05.4)
	The irrigation circuit has not been completely filled	Fill the irrigation circuit via the function FLUSH (see paragraph 05.4)

Problem	Possible cause	Solution
<b>No liquid outflows from the insert during operation</b>	The insert type does not foresee liquid passage	Use an insert type with liquid passage
	The insert is obstructed	Unscrew the insert from the handpiece and free the insert water passage by blowing compressed through it. If the problem persists, replace the insert with a new one
	The handpiece is obstructed	Contact an Authorized Mectron Service Center
	The irrigation level on the screen is adjusted on "0"	Adjust the irrigation level
	The bottle is empty	Fill in the bottle
	The bottle is not correctly installed	Correctly connect the bottle to the device body
	The silicon tubings of the pump are incorrectly installed	Check the connections of the silicon tubings
	The peristaltic pump is over-used	Replace the peristaltic pump (See paragraph 09.1)
<b>Insufficient performance</b>	The insert is not correctly tightened on the handpiece	Unscrew and correctly screw the insert again with the Mectron torque wrench (See paragraph 05.4)
	Insert broken, worn-out, or deformed	Replace the insert with a new one


## 12.3 → REPLACEMENT OF THE FUSES

**⚠ WARNING: Switch the device off.**  
Always turn the device off with the main switch and disconnect it from the electrical power socket before performing the following intervention.

EN

Apply leverage with a flat screwdriver, inserting its tip in the seat of the fuse-holder drawer located under the power supply socket;


1

Pull out the fuse-holder drawer;


**⚠ WARNING: Replace the fuses, complying to the characteristics indicated in Chapter 11 - TECHNICAL DATA**

2

Reinsert the drawer in its housing.

3

## 12.4 → SENDING THE DEVICE TO AN AUTHORISED MECTRON SERVICE CENTRE

In the event that technical assistance is needed on the device, contact one of the Authorized Mectron Service Centers or your Retailer. Do not attempt to repair or modify the device and its accessories.

Clean and sterilize all the parts that need to be sent to an Authorized Mectron Service Center, following the instructions provided in the Chapter 08 - CLEANING AND STERILIZATION. Leave the sterilized parts in the packet, procedure that confirms the sterilization process has been performed.

The cleaning and sterilization requirements are in line with those in force concerning the safeguarding of health and safety on the workplace, as per Law Decree 626/94 and DLgs 81/08 and its subsequent amendments, both laws of the Italian State.

In the case that the customer does not fulfill the requirements, Mectron reserves the right to charge him or her the cleaning and sterilization expenses, or to reject the goods received in unsuitable conditions, returning them to the customer, at his or her expense, for them to be correctly cleaned and sterilized.

The device must be returned properly packaged, accompanied by all the accessories and by a sheet bearing the:

- Data of the owner with telephone number
- Product name
- Serial number and/or lot number
- Reason for goods returned / description of the malfunction
- Photocopy of delivery note or purchase invoice of the device

### **CAUTION: Package**

Pack the device in its original package to prevent damages during transport.

Once the material has been received by the Authorized Mectron Service Center, qualified technical personnel will provide its evaluation to the given circumstance. The repair works will be performed only upon their prior acceptance by the end customer. For further details, contact the closest Authorized Mectron Service Center or your retailer.

Unauthorized repair works may damage the system and make the warranty void, and absolve Mectron of any responsibility for direct or indirect harm to persons or damages to things.

## 13 WARRANTY

EN

All the Mectron devices are subjected to an accurate end inspection that ascertains their full functionality before they are placed on the market.

Mectron provides warranty on the multipiezo purchased new from a retailer or a Mectron importer to cover material or manufacturing defects for:

- 2 YEARS (TWO) on the device from the date of purchase;
- 1 YEAR (ONE) on the handpiece from the date of purchase.

The accessories are not included in the warranty.

During the period in which the warranty is valid, Mectron commits itself to repair (or at its own free choice, replace) those product parts that were to prove defective according to its judgment, free of charge.

Entire replacement of the Mectron products is excluded.

The warranty of the manufacturer and the homologation of the device are not valid in the following cases:

- The device is not used in accordance to the intended use foreseen for it.
- The device is not used in compliance to all the instructions and precautions described in this manual.
- The electrical system in the premises in which the device is used is not compliant to the norms in force and to the relative precautions.

- The assembly operations, extensions, adjustments, updates, and repair works are performed by personnel not authorized by Mectron.

- The environmental conditions for the preservation and storage of the device are not compliant to the precautions indicated in the Chapter 11 - TECHNICAL DATA.

- Use of non-original Mectron inserts, accessories, and spare parts that can compromise correct device operation and cause harm to the patient.

- Accidental breakage due to transport

- Damages due to incorrect use or negligence, or to connection to a power voltage other than the one foreseen

- warranty expiry

**PLEASE NOTE** The warranty is valid only if the warranty sheet annexed to the product has been filled out in all its parts, and sent to our offices or possibly to the retailer or Mectron importer within 20 (twenty) days from the date of purchase, which is the date that appears on the delivery note/purchase invoice issued by the retailer/importer.

To benefit of the warranty service, the customer must return, at his own expense, the device to be repaired to the retailer/Mectron importer from whom he has purchased the product.

See paragraph 12.4 for the details relative to the procedure to send the device to an Authorized Mectron Service Center.

The indications that appear in this publication are not binding and can be modified without fore-notice. The Italian version of this manual is the original document from which its translations have been obtained. In case of any discrepancy, the Italian version will have pertinence.

Texts, images, and graphics of this manual are property of Mectron S.p.A., Carasco, Italy. All rights reserved. The contents cannot be copied, distributed, changed, or made available to third parties without the written approval of Mectron S.p.A.

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