VOCO**Bifluorid 10®**

Fluoride varnish for the treatment of hypersensitivity

Instructions for use

Product description:

Bifluorid 10 is a fast-drying dental suspension that adheres well to dry enamel and dentine. Thus, it is possible to seal the dentine tubules tightly for an extended period of time. **Bifluorid 10** forms a waterproof, protective coat that insulates against thermal and chemical influences.

Indications:

Treatment of all hypersensitivities (dentine hypersensitivity) including:

- hypersensitivities of cervical areas and crown margins
- hypersensitivity after professional cleaning and calculus removal
- hypersensitivity of tooth surfaces after preparation and/or grinding
- hypersensitivity due to worn occlusal surfaces and teeth carrying clasps
- hypersensitivity after injury of the enamel (e.g. fractures, chipping)
- sealing of cavity margins after restoration, especially after the etching technique has been employed
- as dentine protection/lining under amalgam fillings

Method of use:

1. Preparation

Clean and thoroughly air-dry surfaces to be treated.

2. Dispensing

a. Bottle

Shake **Bifluorid 10** bottle well before each use. The solid particles must be well dispersed. The space in the bottle and two metal balls facilitate homogeneous resuspension.

Dispense drops of **Bifluorid 10** onto an applicator or preferably a **Pele Tim** foam pellet. Tightly close **Bifluorid 10** bottle after use.

b. SingleDose

Detach a *SingleDose* application at the perforation and turn the printed side up. Pierce the foil with a **Micro Tim**, enlarge the opening and wet the brush with circling movements.

3. Application

Thinly coat the surfaces. Caution: Layers that are too thick may easily peel. Allow the fluoride varnish to be absorbed for 10 - 20 s, then dry with air. **Bifluorid 10** remains on the applied surfaces for several days if it is correctly applied and appropriate oral hygiene is carried out (**teeth should not be brushed for 12 - 24 h after application**).

For hypersensitive cervical areas, crown margins or similar indications, repeat treatment two to three times at intervals of 7 days as appropriate. Repeated application should only involve the drying of surfaces and reapplication of the varnish.

Information, precautionary measures:

Bifluorid 10 contains ethyl acetate, sodium fluoride, calcium fluoride, cellulose ester and eugenol. **Bifluorid 10** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.

1 g **Bifluorid 10** contains:

50 mg sodium fluoride (equals 23 mg fluoride) and
50 mg calcium fluoride (equals 24 mg fluoride)

Bifluorid 10 should only be used with extreme caution in individuals who may not have full control over their swallowing reflex (e.g. preschool children and the disabled). Systemic fluoride therapy (e.g. fluoride tablets or fluoridated table salt) should be discontinued for several days following application of **Bifluorid 10**.

Incorrect application (see section entitled "Method of use") may result in major amounts of fluoride being introduced into the oral cavity. Swallowing such an amount may lead to nausea, vomiting and / or diarrhoea. Routinely exceeding a total daily amount of 2 mg of fluoride (systemic and local) during odontogenesis can result in disturbances in the mineralisation of the enamel. This disturbance, also referred to as dental fluorosis, manifests itself in the appearance of white spots on the enamel.

The viscosity of **Bifluorid 10** must be low enough to allow the liquid to easily drip from the dropper after being shaken. Add a few drops of the designated **solvent** to thin the preparation if required. The **solvent** should also be used for cleaning the instruments etc.

Storage information:

Protect **Bifluorid 10** from direct sunlight and store at 4 °C - 23 °C. Do not use after the expiry date printed on the package or bottle.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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VC 60 BB1616 E1 0316 99 © by Voco
VOCO GmbH