MM-MTA™

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Activer et mélanger la MM-MTA Caps according to the information in the MM-MTA Caps instruction.

Mixing time for the MM-MTA Caps is 30 seconds.

Attention:

Avoid lag times between the activation, mixing and application processes as the material is in the process of setting which may impair or prevent application of the material. The material must be extruded within 10 seconds after the end of the mixing process. To prevent dehydration during setting, apply MM-MTA™ intraorally immediately after mixing.

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EN Instructions for use

MM-MTA™ – MM-MTA™ Caps Endodontic Repair Cement

MM-MTA™ is an endodontic repair cement. MM-MTA™ powder consists of very fine hydrophilic particles of several mineral oxides. After contact with the liquid, the MM-MTA™ forms a gel that hardens to an impermeable barrier.

MM-MTA™ is delivered in MM-MTA™ Caps. MM-MTA™ Caps are quickly and easily activated and the contents of the capsule are easily ejected with the applicator.

VI. STORAGE

Store MM-MTA™ in a dry place at 10 – 25°C. Do not store below 10°C!

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Capsule mixing (mixing time 30 seconds) is achieved using a high frequency mixer that delivers about 4,300 oscillations/min.

VIII. WARRANTY

S&C Polymer Silicon- und Composite Spezialitäten GmbH warrants this product will be free of material and manufacturing defects. MICRO-MEGA® disclaims all other liabilities and does not implicitly warrant that it is fit for commercialisation or adaptable to a particular purpose. The user shall be liable for determining the suitability of the product for the user's application. In the event of a fault with the product during the warranty period, the sole obligation of S&C Polymer shall be to repair or replace the product of S&C Polymer Silicon- und Composite Spezialitäten GmbH.

IV. INTENDED USES

- Repair of root perforations during root canal treatment.
- Root end filling (retrograde)
- Pulp capping.
- Root end filling (orthograde)

II. CONTRAINDICATIONS

Not known.

III. SIDE EFFECTS

Not known.

IX. LIMITATION OF LIABILITY

Except where prohibited by law, S&C Polymer Silicon- und Composite Spezialitäten GmbH disclaims all liability for any loss or damage arising from this product, whether direct, indirect, specific, incidental or consequential, regardless of the legal grounds, including warranty, contract, negligence or premeditation.

Keep away from children!

For dental use only!

Caution: US Federal law restricts this device to the sale by or on the order of a dentist.

Do not use if package is damaged

Instructions for activating and mixing MM-MTA™ Caps

- Powder
- Liquid
- Plunger

1. MM-MTA™ Cap before activation.

2. For activation of the MM-MTA™ Cap press the plunger on a hard and flat surface until it penetrates the MM-MTA™ Cap.

3. Click once to activate.

Insert the MM-MTA™ Cap into the applicator and click once to activate.

Note: The plunger must level with the rear of the capsule.

4. Insert the MM-MTA™ Cap into a mixer (or an amalgamator), close lid and mix immediately for 30 seconds (about 4300 oscillations / min).

5. Remove the cap from the capsule. If this is not done, the capsule can burst.

6. Insert the MM-MTA™ Cap into the applicator. Pull the lever 2 times (2 clicks) to prime the MM-MTA™ Cap. Extrude the contents of the capsule on a glass plate. Unlock the applicator (push button A) and remove the MM-MTA™ Cap.

The optimal amount of mixed material can only be guaranteed if the applicator is used.

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5.1. Repair of root perforations

Place rubber dam and clean the root canal system using intra-canal instruments and irrigate with NaOCl. Dry the root canal with paper points and isolate the perforation. Fill the apical canal space up to the perforation.

Check the position of MM-MTA™ in the root canal using an X-ray. If the material does not form an even layer, rinse MM-MTA™ out of the canal and repeat the procedure.

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