

RECOMMENDED



OTHER EMS POWDERS

Name		PLUS	CLASSIC COMFORT	CLASSIC	SOFT	PERIO
Ref		DV-082	DV-048	DV-048/ LEM/65	DV-071	DV-070
Flavor		NEUTRAL	LEMON MINT CHERRY NEUTRAL	LEMON	NEUTRAL	NEUTRAL
Composition		ERYTHRITOL	SODIUM BICARBONATE	SODIUM BICARBONATE	GLYCINE	GLYCINE
Particle size		~14µm	~40µm	~65µm	~65µm	~25µm
pH		~7	~8.1	~8.1	~6	~6
Teeth	Stains and early calculus	✓	✓	✓	✓	✗
	Enamel	✓	✓	✓	✓	✓
	Enamel white spots	✓	✗	✗	✗	✓
	Pits and fissures	✓	✗	✗	✗	✗
	Dentine	✓	✗	✗	✗	✓
Soft tissues	Gingiva	✓	✗	✗	✗	✓
	Tongue and palate	✓	✗	✗	✗	✓
Restorations	Hybrid composite	✓	✗	✗	✓	✓
	Glass ionomer	✓	✗	✗	✗	✓
Implants	...and peri-implant tissues	✓	✗	✗	✗	✓
Orthodontics	Brackets and appliances	✓	✗	✗	✗	✓
		✓	✓	✓	✓	✗
		✓	✗	✗	✗	✓
		✓	✗	✗	✗	✓
		✓	✗	✗	✗	✓
		✓	✗	✗	✗	✓

RECOMMENDED



OTHER EMS POWDERS

Name		PLUS	CLASSIC COMFORT	CLASSIC 65 µm	SOFT	PERIO
Ref		DV-082	DV-048/LEM	DV-048/LEM/65	DV-071	DV-070
		✓	✓	✓	✓	✓
		✓	✓	✓	✓	✓
		✓	✓	✓	✓	✓
		✓	✗	✗	✗	✓
		✗	✓	✓	✓	✗
		✓	✓	✓	✓	✓
		✓	✗	✗	✗	✓

► With 2 chemical laboratories in Germany and Switzerland, EMS is the only company producing its own powder. With 35 years of clinical research and testing, many scientific studies and millions of patients treated, EMS recommends the following 2 powders for all clinical cases.