

	No	Type		No	Type		No	Type
AB Dental			Bio3			Conelog		
Narrow Platform	55021	11	Advanced	55067	59	Screw-Line 3.3	55066	58
Standard Platform	55049	41	Progressive	55067	59	Screw-Line 3.8, 4.3	55066	58
						Screw-Line 5.0	55075	67
ADIN			BioHorizons			Cowellmedi		
Touareg S/OS	55034	26	Laser-Lok 3.0 Grey	55034	26	Inno Int Octa	55014	4
Touareg Closefit UNP 2.75	55059	51	Internal Laser-Lok 3.5 Yellow	55049	41	Inno Sub Hex	55037	29
Touareg Closefit NP 3.0	55059	51	Internal Laser-Lok 4.5 Green	55051	43			
Touareg Closefit RP 3.5	55046	38	Internal Laser-Lok 5.7 Blue	55051	43			
Touareg CloseFit WP 4.3, 5.0	55067	59	External 4.0 Platform	55011	1	Dentack		
Swell 3.3	55034	26				Internal Hex	55056	48
			Biomet/3i					
AIRES			External Hex	55011	1	Dental Ratio		
3.3	55049	41	Certain internal Purple	55036	28	Tissue Level RP, WP	55014	4
			Certain internal Purple Platform 3.4	55036	28	Bone Level NC	55061	53
			Certain internal Blue, Yellow, Green	55066	58	Bone Level RC	55062	54
Alfa Gate			T3 Purple	55036	28			
Bioactive, All diam	55034	26	T3 Blue, Yellow, Green	55066	58	Dental Tech		
MAX, All diam	55034	26				Implassic FT3 3.25	55039	31
M+ NP 3.5	55044	36	Bioner			Implassic FT3 3.75	55049	41
M+ RP 4.3 - 5.0	55037	29	HIKELT 3.8-4.7	55034	26			
Porous, All diam	55034	26	Top DM 3.5-5.0	55046	38	Dentalis		
						Regular	55051	43
Alpha Dent			Bionnovation					
Conical connection Narrow Platform	55025	16	Biomorse	55023	14	Dentegris		
Conical connection Wide Platform	55057	49				3.3	55034	26
Hexagon connection Narrow Platform	55039	31	Biotech Dental			3.75	55057	49
Hexagon connection Regular Platform	55057	49	Kontakt, Kontakt-S 3.6 - 4.8	55045	37	4.5	55049	41
						3.8/4.5	55049	41
Alpha-Bio			BlueSkyBio			5.5	55020	10
Internal Hex	55034	26	BIO Max NP Platform, All	55046	38			
Conical Hex 3.2	55065	57	BIO Max RP Platform 4.3-5.0	55037	29	Dentis		
Conical Hex 3.75-5.0	55051	43	BIO Int Hex 3.5 Platform 3.7-4.1	55049	41	Cleanlant S-Clean 3.7 - 4.8	55020	10
			BIO Int Hex 3.5 Platform 4.7	55051	43	S-Clean OneQ Narrow 3.3	55064	56
Anker			BIO Int Hex 4.5 Platform 4.7-7.0	55051	43			
SB III 3.5	55023	14	BIO One-stage Narrow, Regular, Wide	55014	4	Dentium		
SB III 4.0-5.0	55020	10	BIO One-stage Short Regular & Wide	55014	4	Superline	55037	29
			BIO Quattro Narrow 3.3	55061	53	NR Line	55066	58
Ankylos			BIO Quattro Regular 4.1 - 5.6	55062	54	Simple Line II GH 1,6	55014	4
C/X	55025	16	BIO Conus 12, 3.0 Platform 3.0	55064	56			
			BIO Conus 12, 3.5/4.0 Platform 3.5-4.0	55046	38	DEVA		
Anthogyr			BIO Conus 12, 4.5 x 6mm	55046	38	Bone Level 3.5	55064	56
Axiom Bone Level 3.4 - 5.2	55065	57	BIO Conus 12, 4.5/5.0 Platform 4.5 - 5.0	55017	7	Bone Level 4.5	55051	43
Axiom Tri-lobe	55065	57	BIO Conus 12, 5.0 x 6mm	55046	38			
Axiom Tissue Level Narrow, Regular	55050	42	BIO Trilobe 3.5	55016	6	DIO		
			BIO Trilobe 4.3-6.0	55017	7	SM Submerged 3.8	55049	41
Argon						SM Submerged 5.0	55020	10
K3pro Rapid 3.0 - 3.5	55026	17	Bredent			UF II Narrow	55064	56
K3pro Rapid Platf 3mm	55019	9	SKY	55025	16	UF II Regular	55020	10
K3pro Rapid Platf 3mm Short	55044	36	copaSKY	55027	19			
			mini2SKY	55033	25	Ditron Dental		
Astra Tech			fast & fixed Abutment	55083	75	Ultimate 3.75-6.0	55087	79
Osseospeed EV 3.0	55058	50				MPI 3.3-6.0	55087	79
Osseospeed EV 3.6	55046	38	BTI					
Osseospeed EV 4.2	55057	49	Interna Universal	55051	43	DSP Biomedical		
Osseospeed EV 4.8	55055	47	Interna Wide	55051	43	External Hexagon (HE) Standard	55014	4
Osseospeed EV 5.4	55048	40				Hexagonal Internal Standard (HiS)	55034	26
Osseospeed TX 3.0 S	55064	56	BTK			Hexagonal Internal Morse (HiM)	55046	38
Osseospeed TX 3.5-4.0 S	55046	38	External Hex EN	55038	30	Hexagonal Internal Large (HiL)	55037	29
Osseospeed TX 4.5-5.0 S	55017	7	External Hex ER	55011	1	Indexed Morse Taper (CMI)	55025	16
			External Hex EW	55011	1			
Avinent			Internal Hex IR	55049	41	Easy Implant		
External, Coral & Ocean 3.5	55038	30	Internal Hex IM	55051	43	Ocean	55046	38
External, Coral & Ocean 4.1	55011	1	Internal Hex IW	55051	43	Lilas	55037	29
External, Coral 5.1	55013	3	Hexagonal Morse Taper KR	55081	73			
Internal, Coral & Ocean 3.5, 4.1	55051	43	Hexagonal Morse Taper KW	55056	48	Edierre		
Conical, Ocean 3.5	55066	58	Octagonal Morse Taper SR	55014	4	3-lobe	55021	11
Conical, Ocean 4.1	55004	18				4-lobe	55025	16
			C-Tech					
B&B Implant			EL	55067	59	Euroteknika		
Duravit Slim	55021	11				Internal Platform 3.0	55064	56
Duravit 3P	55035	27	Camlog			Internal Platform 3.5	55081	73
Duravit Evolution	55035	27	Screw-Line, Root-Line 3.3	55036	28	Naturactis/Naturall+ 3.0	55064	56
			Screw-Line, Root-Line 3.8, 4.3	55031	23	Naturactis 3.5-5.0	55081	73
BEGO			Screw-Line, Root-Line 5.0, 6.0	55032	24	Uneva	55011	1
Semados Mini	55012	2	iSy	55065	57			
Semados S / RI 3.25-5.5	55034	26				galimplant		
Semados RS / RSX 3.0	55082	74	Ceralog			iPx 4012	55046	38
Semados RS / RSX 3.75-5.5	55034	26	Hexalobe ceramic	55076	68			
Semados SC / SCX 3.25-5.5	55034	26						

	No	Type		No	Type		No	Type
GC Tech			Implant One			Leader		
Aadva Narrow	55081	73	100 Series Std thread	55029	21	Internal Hex, 3.3 - 5.5 mm	55034	26
Aadva Regular	55066	58	200 Series Std thread	55029	21	Cono-In	55067	59
			300 Series Std thread (blue)	55029	21			
Geass			400 Series Std thread (purple)	55035	27	Maxillent		
Way Extra	55081	73	500 series Std thread (green)	55035	27	iRaise	55049	41
			300 Series Wide thread (blue)	55029	21			
Global D			400 Series Wide thread (purple)	55035	27	Medentika		
In-Kone 4.0-5.0	55024	15	500 Series Wide thread (green)	55035	27	Microcone NI	55021	11
TwinKon 3.5-4.5	55122	88				Microcone RI	55081	73
			Implant Swiss			Quattrocone RI	55081	73
gmi			Tissue level	55014	4	Quattrocone 30 AI	55081	73
avantgard RP 3.75-4.25	55066	58	Bone level Regular	55037	29	Procone 3.3	55119	85
avantgard WP 4.75-5.75	55075	67	Bone level 3.3	55065	57	Procone 3.8	55031	23
frontier RP,WP	55051	43				Procone 4.3	55031	23
phoenix NP 3.3	55038	30	Implantvision			Procone 5.0	55032	24
phoenix Std 3.3 - Wide 5.0	55011	1	Internal Hex 3.7-5.7	55034	26			
Monolith 3.0	55043	35				Medentis		
			Innobiosurg			ICX Premium 3.75 - 4.8	55062	54
GP Implants			MagiCore Regular	55018	8	ICX Plus 3.45	55039	31
SSI Spiral Shape implant	55051	43				ICX Narrow	55058	50
			Intoss Anchor					
gt medical			Internal Hex 3.7-5.7	55088	80	Medical Instinct		
HXE	55011	1				Hex 3.0	55071	63
HXE RP	55012	2	Intra-lock			Hex 3.4, 4.0, 5.0	55052	44
HXE WP	55013	3	Intrahex 3.4 - 4.5	55051	43	Cone 3.0, 3.4, 4.0	55071	63
HXI 4.1	55037	29	Gold & Blue 3.4 - 6.0	55051	43	Cone 5.0	55072	64
002 HXI 4.7	55049	41	InDex 3.4 - 3.75	55015	5			
OCT 4.3	55014	4	InDex 4.0 - 4.3	55020	10	Medigma Biomedical		
			FlatTop 3.75	55012	2	Mars 3,5-6,0	55057	49
Hahn						FixTite 3,5-6	55057	49
Tapered 3.0	55058	50	IRES			Fix-a-Dent 3,75-6	55057	49
Tapered 3.5	55046	38	External Hex Narrow 3.3	55012	2			
Tapered 4.3	55046	38	External Hex Regular 3.7 - 5.2	55011	1	Medimecca		
Tapered 5.0	55037	29	Internal Hex Narrow 3.3	55087	79	PT SLA Chaorum	55047	39
Tapered 7.0	55048	40	Internal Hex Regular 3.7 - 5.2	55049	41			
			Internal Octagon Regular 3.7 - 5.2	55014	4	Megagen		
Hi-Tec						Anyridge	55051	43
Tapered Self Thread NTI 3.3, SLT 3,75	55034	26	J Dental Care			Anyone Internal	55037	29
Tapered Self Thread GTI 4.2, WTI 5.0	55051	43	JDEvolution 3.7-6.0	55067	59	Mini	55064	56
Spark TPI	55034	26	JDEvolution Plus+ 3.7-6.0	55067	59	MUA	55080	72
Logic Plus	55034	26	JD Evolution S 3.25	55044	36	Octa	55121	87
Logic LGI 3.0-3.5	55046	38	JD ICON ultra, S 2.75	55064	56			
Logic LGI 4.3-5.0	55037	29	Conical Abutment	55080	72	Microdent		
						Universal	55011	1
Hiossen			Keystone			Microdent 2.8 - 3.25	55068	60
HS System	55028	20	Prima / Connex	55049	41	Microdent 3.3 - 5.0	55011	1
ET System	55037	29	Prima / Connex	55067	59	Microdent 5.5 - 8.0	55013	3
ET 3.2	55047	39	Prima / Connex	55041	33	Genius 3.0	55025	16
ET mini	55015	5	Prima / Genesis Narrow	55049	41	Genius 3.5 - 5.0	55051	43
			Prima / Genesis Regular	55067	59	Ektos 3.2	55023	14
IDI Evolution			Prima / Genesis Wide	55041	33	Ektos 3.7 - 8.0	55051	43
Standard implant 3.75-4.0	55004	18	Restore Narrow	55079	71	Trylogic 3.5 - 4.2	55051	43
B1One High	55080	72	Restore Regular	55011	1	Trylogic 5.0	55075	67
B1One Low	55033	25	Restore Wide	55013	3			
			Klockner			MIS		
IDI Implants			VEGA Micro 3.0	55074	66	Conical connection Narrow (C1, V3)	55081	73
BIO 3.5-5.0	55037	29	VEGA Narrow 3.5	55065	57	Conical connection Standard (C1, V3)	55087	79
ID Cam ST 4.2	55010	13	VEGA Regular 4.0, 4.5	55034	26	Conical connection Wide (C1)	55041	33
ID All	55010	13	KL Narrow 3.3	55012	2	Seven and M4 Narrow	55027	19
ID3 3.1	55021	11	KL Regular 3.7, 4.2	55011	1	Seven and M4 Standard	55049	41
			KL Wide 4.7	55013	3	Seven Wide 5.0	55051	43
ids						Connect abutment NP, SP, WP	55115	81
Reflect Aspire	55081	73						
Reflect Certus	55019	9						
Reflect Tapered	55034	26						
Reflect Recover	55020	10						
Reflect Rapid	55020	10						
			Konmet			myplant GmbH		
Implance			Classic 4.0, 4.8	55014	4	Myplant two 3,5-4,5	55077	69
Tissue level	55014	4	Cylindrical 3.5-4.0	55065	57			
Bone level Regular	55037	29	Cylindrical 4.5-5.5	55057	49	Neo Biotech		
Bone level 3.3	55065	57	Short Cylindrical 5.0-6.0	55078	70	IS II 4.0	55020	10
			Conical 4.4/3,4	55065	57	IS III 4.0	55020	10
			Conical 4.8/3,8, 5,5/4,5	55057	49	Neodent		
Implant Direct						Drive CM	55025	16
Legacy 1/2/3/4 3.0/3.2	55012	2	Kyocera			Titamax WS	55056	48
Legacy 1/2/3/4 3.5/3.7 - 5.7/7.0	55051	43	FINESIA Bone Level 3.2-3.4 (NP)	55058	50	GM	55066	58
RePlant 3.5	55022	12	FINESIA Bone Level 3.7-4.2 (RP)	55034	26			
RePlant 4.3 - 5.0	55010	13	FINESIA Tissue Level 3.7-4.7 (RP)	55014	4	Neoss		
RePlant 6.0	55032	24	FINESIA Tissue Level 4.7-5.2 (WP)	55014	4	SP (Ø 3.5 implants and wider)	55004	18
ReActive 3.5 Platform	55022	12				NP (Ø 3.25)	55027	19
ReActive 4.3 - 5.0 Platform	55010	13	Lasak			Access Abutment	55060	52
Interactive 3.0 Platform	55046	38	BioniQ QR	55066	58			
Interactive 3.4 Platform	55037	29	BioniQ QN	55021	11			
Screwindirect 5.0 PF 3.0 - 5.7	55070	62						

	No	Type		No	Type		No	Type
Nobel Biocare			Paltop			SIN		
Brånemark System Ext hex RP	55011	1	Standard 3.75-5.0	55034	26	Strong SW CM	55023	14
Brånemark System Ext hex NP	55012	2	Narrow	55027	19	Strong SW Int Hex	55035	27
Brånemark System Ext hex WP	55013	3	PCA Conical	55066	58	Unitite CM 2.9	55084	76
Tri-lobe Connection NP	55022	12				Unitite CM 3.5	55077	69
Tri-lobe Connection RP, WP	55010	13	Phibo			Unitite CM Compact 4.0	55085	77
Tri-lobe Connection 6.0	55032	24	Aurea EVO NP	55027	19	Tryon 3.25 - 5.0	55011	1
Conical Connection 3.0	55058	50	Aurea EVO RP, WP	55051	43			
Conical Connection Narrow	55046	38	TSA S3	55074	66	Southern Implants		
Conical Connection Regular	55037	29	TSA S4, S5	55034	26	Tri-Nex 4.3-5.0	55010	13
Conical Connection Wide 5.5	55048	40	TSH S3, S4, S5	55011	1	Tri-Nex 3.5	55022	12
MUA NP & RP	55080	72				External Hex	55011	1
On1 NP, RP, WP	55042	34	P-I Brånemark			MAX Ext Hex	55011	1
			HEX 3.3	55115	81	Co-Axis Ext Hex	55011	1
			HEX 3.75 -5.0	55117	83	Piccolo Ext Hex	55050	42
Noris Medical			AMP 3.3	55039	31	IT-Connection	55014	4
Tuff 3.3 - 6.0	55034	26	AMP 3.75 - 4.8	55118	84	Deep Conical 3.0	55064	56
			MT 3.3 - 4.8	55046	38	Deep Conical 3.5-4.0	55046	38
Novem			Conical Abutment HEX, AMP, MT	55080	72	Deep Conical Tapered 5.0	55017	7
3.4	55039	31						
3.8	55016	6						
4.2, 5.0	55057	49	Prodent Italia			Spiral Tech		
			Twinner & Prime, Conomet, Collar: Orange	55027	19	ESi, Hex connection	55051	43
			Twinner & Prime, Conomet, Collar: Fuchsia	55051	43	Ultimate, Hex connection	55051	43
NucleOss			Twinner & Prime, Conomet, Collar: Green	55067	59	ESi Conical 3.5	55044	36
Tpure Bone Level	55004	18	Twinner & Prime, Conomet, Collar: Yellow	55020	10	ESi Conical 4.3, 5.0	55037	29
T6 Bonel Level	55074	66	Prime, Conomet, Collar: Blue	55020	10			
						Straumann		
OCO Biomedical						Tissue level, Regular & Wide Neck	55014	4
Engage 3.25	55027	19	Qualia			Tissue level, Narrow Neck Crossfit (NNC)	55076	68
Engage 4.0-6.0	55051	43	NeoActive	55046	38	Bone level / Bone Level Tapered Narrow	55061	53
ERI 3.25	55069	61				Bone level / Bone Level Tapered Regular	55062	54
ERI 4.0-5.0	55014	4	Quest			Bone level / Bone Level Tapered SC 2.9	55064	56
TSI 3.25	55068	60	Hylock NP	55030	22	Bone Level X (BLX) all platform	55027	19
TSI 4.0-5.0	55011	1	Hylock RP	55066	58	Screw retained Abutment NP, RP	55080	72
OsseoFuse			Resista			Surgikor		
Hexaplus/ HexaPlus-S Sinus 3.75-7.5	55081	73	WIS 3.0-5.0	55049	41	Fixation 3.0	55058	50
ConicalPlus/ ConicalPlus-S Sinus 3.75-5.25	55015	5				Fixation 3.5	55066	58
			Rex Implants			Fixation 4.3	55037	29
			External hex TL 1.8	55060	52	Solution 3.5, 4.0	55066	58
Osstem			Internal hex TL 2.9	55034	26	Solution 4.5, 6.0	55037	29
US Regular	55011	1				Immediate, all	55049	41
US Mini	55012	2	Ritter			Versatile, all	55049	41
US Wide	55013	3	Multi Unit	55033	25			
US Wide PS	55011	1	Spiral SB/LA Narrow PF 2.9	55081	73	Sweden & Martina		
SS Regular, Wide	55028	20	Spiral SB/LA Standard PF 3.75	55049	41	Syra, Syra SP 3.9-5.0	55011	1
TS Regular, Reg Short, Diam 6 & 7	55037	29				Outlink ² 3.75, 4.10, 5.0	55011	1
TS Mini 3.5	55015	5	RS Implant			Outlink ² 3.3	55079	71
TS Mini 3.0	55047	39	3.5 - 6.0	55020	10	Shelta 3.8-5.0	55025	16
						Premium Straight & SP 3.3-3.8	55025	16
Osteo Plus			Schütz Dental			Kohno Straight & SP 3.3-3.8	55025	16
Iridium PF 3.8 & 4.2	55079	71	Impla cylindrical	55052	44	Premium Straight & SP 4.25-6.0	55004	18
She 3.45, 3.75	55011	1	IMPLA micro Retention Cone 4.2	55052	44	Kohno Straight & SP 4.25-6.0	55004	18
Shi 3.75, 4.2	55049	41				Prama 3.8 - 5.0	55079	71
Shi 5.0	55051	43	SGS Implant			DAT Connection	55035	27
			P1, P7 Platform	55049	41			
OT-Medical			P1D, P7D	55049	41	Tag Dental		
OT-F2 3.4, 3.8	55053	45	P7N	55059	51	Axis	55051	43
OT-F2 4.1	55073	65				massif	55051	43
			SIC invent					
Overmed			MUA	55120	86	Tatum		
Easyline Easydip	55020	10	SiCace 3.4 - 5.0	55052	44	Integrity 3.7 - 5.0	55015	5
ETLINE	55034	26	SiCmax, SiCtapered 3.7 - 5.2	55052	44	Integrity 6.0 - 8.0	55040	32
			SiCvantage® max, SiCvantage tapered 3.0	55026	17	T implant	55063	55
Oxy Implant			SiCvantage® max, SiCvantage tapered 3.7	55116	82			
Internal	55020	10	SiCvantage® max, SiCvantage tapered 4.2-5.2	55054	46	TBR		
Kone Conical Narrow	55021	11				Infinity PF 3.5, PF 4.0	55004	18
Kone Conical Regular	55051	43	Simply Implants			M Implant	55044	36
PSK Line Regular	55051	43	Internal Hex	55051	43			
External Mini	55079	71				Thommen		
External Regular	55011	1				SPI 3.5	55021	11
External Maxi	55011	1				SPI 4.0	55065	57
						SPI 4.5-6.0	55066	58

	No	Type
ticare		
Osseus Mini, Standard	55011	1
Osseus Maxi	55013	3
Osseus Maxi ps	55011	1
Inhex Mini	55021	11
Inhex Standard	55046	38
Inhex Maxi	55017	7
Osseus Quattro Standard	55011	1
Inhex Quattro Standard	55046	38
TRATE		
ROOTFORM	55065	57
TRI Dental implants		
Tissue level Octa	55014	4
Tri-Vent, all	55034	26
Narrow	55025	16
Trinon		
Q2-Implant	55011	1
QK-Implant	55014	4
QZA 3.5 Platform	55049	41
QZA 4.5, 5.7 Platform	55051	43
URIS		
Omni Narrow 3.0-3.5	55115	81
Omni Regular 4.0-6.5	55020	10
Vulkan		
Internal Hex	55034	26
External Hex NP	55012	2
External Hex RP	55011	1
Conical con MiP	55058	50
Conical con NP	55046	38
Winsix		
Internal Hex	55020	10
Xgate		
Conical 3.3 - 5.0	55034	26
Xive		
S-Plus 3.0	55030	22
S-Plus 3.4	55061	53
S-Plus 3.8-5.5	55062	54
TG Plus 3.4-4.5	55018	8
Zeramex		
XT Implant 4.2, 5.5	55086	78
Ziacom		
Conical connection RP	55046	38
Conical connection WP	55075	67
External hexagonal connection NP	55079	71
External hexagonal connection RP	55117	83
External hexagonal connection WP	55014	4
Internal hexagonal connection NP	55027	19
Internal hexagonal connection RP	55057	49
Internal hexagonal connection WP	55057	49
Zimmer		
Screw-Vent, Tapered Screw-Vent 3.7-4.1	55049	41
Screw-Vent, Tapered Screw-Vent 4.7-5.7	55051	43
Tapered Swiss Platf 3.8	55070	62
Tapered Swiss Platf 4.8	55014	4
Spline Twist Platf 4 & 5	55038	30
Eztetic NP 3.1	55065	57
Zuga		
All	55027	19